



DAG

SURGICAL CORPORATION



Liposuction

CANNULAS & ACCESSORIES

✉ info@dagsurgical.co

☎ +92 3498306020

📍 DAG SURGICAL CORPORATION
Haji Pura, Church Road,
Sialkot-51310 Pakistan.

www.dagsurgical.co

Welcome to DAG SURGICAL CORPORATION

Introduction DAG SURGICAL CORPORATION sustained an enviable reputation since 1990 being the manufacturer of precision surgical instruments covering each type of surgical specialty. We have distinguished ourselves by providing the high-quality surgical instruments and gained the loyalty, support and innovative ideas from our precious customers worldwide.

Vision: DAG SURGICAL CORPORATION emphasis on an unparalleled experience being a trusted individual for business partners, health care professionals and the humanity.

Mission:
Being reliable manufacturer supply for all the surgical requirements by delivering quality at an affordable price with great services and experience.

About Us: DAG SURGICAL CORPORATION was established to serve the medical industry with quality surgical instruments under the experience and expertise of many craftsmen. Where traditional manufacturing skills are combined with modern manufacturing techniques, resulting in precious quality instruments to meet and enhance customer's expectations level. With our flexibility and ability to accomplish an individual's needs we earned the confidence of our clients worldwide, which enable us to sustain long-term business relationships with them, based on principles of trust, collaboration, open communication, professional delivery and innovation.

Production:
We have in-house production to work with forgings and dies. And being a manufacturer, we not only manufacture standard products, but also are capable to serve with a range of services and technologies that enable us to work on custom forgings, drawings and prototypes under assured accuracy.

Quality Certifications: We are certified:
1 - ISO-9001:2015
2 - ISO-13485:2016
3 – Medical Device Directive 93/42/EEC
4 – FDA - USA

Markets:
Our main market of exports includes:


- Europe
- Germany
- USA
- Turkey

Equal Employment Opportunity:
We are committed to equal employment opportunity regardless of race, color, ancestry, religion, sex, national origin, sexual orientation, age, citizenship, marital status, disability, gender identity or Veteran status.

Table of Contents

- Cannula Hubs/
Connections -
MICROAIRE-TYPE
POWER CANNULA- -
Infiltration Cannula -
General Suction
Cannula - Harvester
Cannula - Injector
Cannula - Speciality
Cannula - Accessories

Cannula Hubs/Connections



Luer Lock Hub/Connector:

The universal Luer Lock hub enables cannula to attach directly with any Luer Lock Syringe.

Luer Lock Stainless Steel Hub/Connector:

The universal Luer Lock hub enables cannula to attach directly with any Luer Lock Syringe.

Super Luer Lock Hub/Connector:

A modification to general Luer Lock hub which is more powerful, more stable, Leak proof, free from instability and contamination.

Threaded Hub/Connector:


Threaded connector hub style is used mostly for machine assisted Liposuction cannulas, Dag Surgical corp's unique Light weight threaded cannula can also be used with ordinary syringe. All you need is Dag Surgical corp's connector from threaded to Luer style.

60cc Tommey Hub/Connector:

An ideal hub style for surgeons why love to handle large fat volumes. This hub style connect cannula to a 60cc Tommey to syringe.

Fixed Handle Cannulas:

DAG SURGICAL CORPORATION super power grip handle fixed with the cannula of your own choice! The unique design of handle offer powerful grip harvesting large volumes. You will feel amazingly comfortable ever after a long operational time. The thumb grip provide a unique opportunity to control the thrust.



MICROAIRE-TYPE POWER CANNULA



MICROAIRE-TYPE POWER CANNULA



DAG SURGICAL CORPORATION has developed a superior, reusable cannula, fabricated from the highest-quality materials, using tight tolerances to minimize vibration. Now available in a full line of sizes and patterns to meet all your liposuction needs.

- High quality materials and craftsmanship • Cannula tubing connector designed for a non-slip fit for both aspiration and infiltration tubing • Reusable / Autoclavable • Validated cleaning and sterilization instructions

CANNULA PATTERNS



ACCELERATOR

DAG-III-326 PWR	3MM X 26CM STRAIGHT
DAG-III-426 PWR	4MM X 26CM STRAIGHT
DAG-III-432 PWR	4MM X 32CM STRAIGHT
DAG-III-436 PWR	4MM X 36CM STRAIGHT
DAG-III-526 PWR	5MM X 26CM STRAIGHT
DAG-III-532 PWR	5MM X 32CM STRAIGHT



BASKET

DAG-BAS-326 PWR	3MM X 26CM STRAIGHT
DAG-BAS-332 PWR	3MM X 32CM STRAIGHT
DAG-BAS-426 ANG PWR	4MM X 26CM ANGLED
DAG-BAS-426 PWR	4MM X 26CM STRAIGHT
DAG-BAS-432 ANG PWR	4MM X 32CM ANGLED
DAG-BAS-432 PWR	4MM X 32CM STRAIGHT
DAG-BAS-440 ANG PWR	4MM X 40CM ANGLED
DAG-BAS-440 PWR	4MM X 40CM STRAIGHT
DAG-BAS-526 ANG PWR	5MM X 26CM ANGLED
DAG-BAS-526 PWR	5MM X 26CM STRAIGHT
DAG-BAS-532 ANG PWR	5MM X 32CM ANGLED
DAG-BAS-532 PWR	5MM X 32CM STRAIGHT
DAG-BAS-536 ANG PWR	5MM X 36CM ANGLED
DAG-BAS-536 PWR	5MM X 36CM STRAIGHT



DOUBLE BASKET

DAG-DBAS-426 PWR	4MM X 26CM STRAIGHT
DAG-DBAS-432 ANG PWR	4MM X 32CM ANGLED
DAG-DBAS-432 PWR	4MM X 32CM STRAIGHT
DAG-DBAS-440 ANG PWR	4MM X 40CM ANGLED
DAG-DBAS-440 PWR	4MM X 40CM STRAIGHT
DAG-DBAS-526 PWR	5MM X 26CM STRAIGHT
DAG-DBAS-532 PWR	5MM X 32CM STRAIGHT



CANDY CANE

DAG-CC-326 PWR	3MM X 26CM STRAIGHT
DAG-CC-426 PWR	4MM X 26CM STRAIGHT
DAG-CC-432 PWR	4MM X 32CM STRAIGHT
DAG-CC-436 PWR	4MM X 36CM STRAIGHT
DAG-CC-526 PWR	5MM X 26CM STRAIGHT
DAG-CC-532 PWR	5MM X 32CM STRAIGHT



BLACK HOLES

DAG-CC-326 PWR	5MM X 15CM STRAIGHT
DAG-CC-426 PWR	5MM X 20CM STRAIGHT
DAG-CC-432 PWR	5MM X 25CM STRAIGHT
DAG-CC-436 PWR	5MM X 30CM STRAIGHT
DAG-CC-526 PWR	5MM X 35CM STRAIGHT



MERCEDES

DAG-MER-326 PWR	3MM X 26CM STRAIGHT
DAG-MER-426 ANG PWR	4MM X 26CM ANGLED
DAG-MER-426 PWR	4MM X 26CM STRAIGHT
DAG-MER-432 ANG PWR	4MM X 32CM ANGLED
DAG-MER-432 PWR	4MM X 32CM STRAIGHT
DAG-MER-436 ANG PWR	4MM X 36CM ANGLED
DAG-MER-436 PWR	4MM X 36CM STRAIGHT
DAG-MER-526 ANG PWR	5MM X 26CM ANGLED
DAG-MER-526 PWR	5MM X 26CM STRAIGHT
DAG-MER-532 ANG PWR	5MM X 32CM ANGLED
DAG-MER-532 PWR	5MM X 32CM STRAIGHT
DAG-MER-536 ANG PWR	5MM X 36CM ANGLED
DAG-MER-536 PWR	5MM X 36CM STRAIGHT



DOUBLE MERCEDES

DAG-DMER-326 PWR	3MM X 26CM STRAIGHT
DAG-DMER-426 PWR	4MM X 26CM STRAIGHT
DAG-DMER-432 PWR	4MM X 32CM STRAIGHT
DAG-DMER-436 PWR	4MM X 36CM STRAIGHT
DAG-DMER-526 PWR	5MM X 26CM STRAIGHT
DAG-DMER-532 PWR	5MM X 32CM STRAIGHT
DAG-DMER-536 PWR	5MM X 36CM STRAIGHT



HELIX 3 PORT

DAG-HLX-326 PWR	3MM X 26CM STRAIGHT
DAG-HLX-332 PWR	3MM X 32CM STRAIGHT
DAG-HLX-426 PWR	4MM X 26CM STRAIGHT
DAG-HLX-432 PWR	4MM X 32CM STRAIGHT
DAG-HLX-436 PWR	4MM X 36CM STRAIGHT
DAG-HLX-526 PWR	5MM X 26CM STRAIGHT
DAG-HLX-532 PWR	5MM X 32CM STRAIGHT
DAG-HLX-536 PWR	5MM X 36CM STRAIGHT



LARGE VOLUME HARVESTER

DAG-LVH-326 PWR	3MM X 26CM STRAIGHT
DAG-LVH-332 PWR	3MM X 32CM STRAIGHT
DAG-LVH-426 PWR	4MM X 26CM STRAIGHT
DAG-LVH-432 PWR	4MM X 32CM STRAIGHT
DAG-LVH-440 PWR	4MM X 40CM STRAIGHT
DAG-LVH-526 PWR	5MM X 26CM STRAIGHT
DAG-LVH-532 PWR	5MM X 32CM STRAIGHT



SPATULA

DAG-SPA-315 PWR	3MM X 15CM STRAIGHT
DAG-SPA-415 PWR	4MM X 15CM STRAIGHT



BLACK HOLES

DAG-CC-326 PWR	3MM X 15CM STRAIGHT
DAG-CC-426 PWR	3MM X 20CM STRAIGHT
DAG-CC-432 PWR	3MM X 25CM STRAIGHT
DAG-CC-436 PWR	3MM X 30CM STRAIGHT
DAG-CC-526 PWR	3MM X 35CM STRAIGHT



LUER LOCK CANNULA

(Stainless Steel)



CANNULAS IN DIFFERENT SHAPES



Infiltration Cannulas

DAG SURGICAL CORPORATION Infiltration Cannulas are designed for safe and efficient delivery of tumescent solution. The spiral micro holes on cannula ensure and even distribution of local anesthetic which results in a longer anesthetic effects this less painful and more pleasant liposuction experience for patient.




**CANNULA
IN BLACK**

1213



Infiltration Cannulas



Tumescant Infiltrator

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60ccTommeyhub

DAG-101.01	Ø 2.5 mm X 15 cm
DAG-101.02	Ø 2.5 mm X 20 cm
DAG-101.03	Ø 3.0 mm X 20 cm
DAG-101.04	Ø 3.0 mm X 25 cm

Treaded Hub

DAG-102.01	Ø 2.0 mm X 15 cm
DAG-102.02	Ø 3.0 mm X 15 cm
DAG-102.03	Ø 3.0 mm X 20 cm
DAG-102.04	Ø 3.0 mm X 30 cm
DAG-102.05	Ø 4.0 mm X 20 cm
DAG-102.06	Ø 4.0 mm X 30 cm

Ultra Luer lock

DAG-103.01	Ø 1.0 mm X 5 cm
DAG-103.02	Ø 1.0 mm X 7 cm
DAG-103.03	Ø 1.5 mm X 7 cm
DAG-103.04	Ø 1.5 mm X 9 cm

Luer lock

DAG-104.01	Ø 1.0 mm X 4 cm
DAG-104.02	Ø 1.0 mm X 5 cm
DAG-104.03	Ø 1.5 mm X 7 cm
DAG-104.04	Ø 1.5 mm X 9 cm
DAG-104.05	Ø 1.5 mm X 10 cm
DAG-104.06	Ø 1.5 mm X 12 cm
DAG-104.07	Ø 2.0 mm X 10 cm
DAG-104.08	Ø 2.0 mm X 12 cm
DAG-104.09	Ø 2.0 mm X 15 cm
DAG-104.10	Ø 2.0 mm X 20 cm
DAG-104.11	Ø 2.5 mm X 15 cm
DAG-104.12	Ø 2.5 mm X 20 cm
DAG-104.13	Ø 3.0 mm X 15 cm
DAG-104.14	Ø 3.0 mm X 20 cm



Luer lock (Stainless Steel)

DAG-105.01	Ø 1.0 mm X 4 cm
DAG-105.02	Ø 1.0 mm X 5 cm
DAG-105.03	Ø 1.0 mm X 7 cm
DAG-105.04	Ø 1.5 mm X 7 cm
DAG-105.05	Ø 1.5 mm X 9 cm
DAG-105.06	Ø 1.5 mm X 10 cm
DAG-105.07	Ø 1.5 mm X 12 cm
DAG-105.08	Ø 2.0 mm X 10 cm
DAG-105.09	Ø 2.0 mm X 12 cm
DAG-105.10	Ø 2.0 mm X 15 cm
DAG-105.11	Ø 2.0 mm X 20 cm
DAG-105.12	Ø 2.5 mm X 15 cm
DAG-105.13	Ø 2.5 mm X 20 cm
DAG-105.14	Ø 3.0 mm X 15 cm
DAG-105.15	Ø 3.0 mm X 20 cm

Power Cannula

DAG-106.01	Ø 2.0 mm X 15 cm
DAG-106.02	Ø 2.0 mm X 20 cm
DAG-106.03	Ø 2.0 mm X 25 cm
DAG-106.04	Ø 2.0 mm X 30 cm
DAG-106.05	Ø 2.0 mm X 35 cm
DAG-106.06	Ø 3.0 mm X 15 cm
DAG-106.07	Ø 3.0 mm X 20 cm
DAG-106.08	Ø 3.0 mm X 25 cm
DAG-106.09	Ø 3.0 mm X 30 cm
DAG-106.10	Ø 3.0 mm X 35 cm
DAG-106.11	Ø 4.0 mm X 15 cm
DAG-106.12	Ø 4.0 mm X 20 cm
DAG-106.13	Ø 4.0 mm X 25 cm
DAG-106.14	Ø 4.0 mm X 30 cm
DAG-106.15	Ø 4.0 mm X 35 cm
DAG-106.16	Ø 5.0 mm X 15 cm
DAG-106.17	Ø 5.0 mm X 20 cm
DAG-106.18	Ø 5.0 mm X 25 cm
DAG-106.19	Ø 5.0 mm X 30 cm
DAG-106.20	Ø 5.0 mm X 35 cm
DAG-106.21	Ø 6.0 mm X 15 cm
DAG-106.22	Ø 6.0 mm X 20 cm
DAG-106.23	Ø 6.0 mm X 25 cm
DAG-106.24	Ø 6.0 mm X 30 cm
DAG-106.25	Ø 6.0 mm X 35 cm

Infiltration Cannulas



Tumescant Bayonet Infiltrator

The bayonet shaped handles allow an easier manipulation

60ccTommeyhub

DAG-101.11	Ø 2.0 mm X 15 cm
DAG-101.12	Ø 2.5 mm X 20 cm
DAG-101.13	Ø 3.0 mm X 20 cm
DAG-101.14	Ø 3.0 mm X 25 cm

Treaded Hub

DAG-102.11	Ø 2.0mm X 15 cm
DAG-102.12	Ø 3.0mm X 15 cm
DAG-102.13	Ø 3.0mm X 20 cm
DAG-102.14	Ø 3.0mm X 30 cm
DAG-102.15	Ø 4.0mm X 20 cm
DAG-102.16	Ø 4.0mm X 30 cm

Ultra Luer lock

DAG-103.11	Ø 1.0 mm X 5 cm
DAG-103.12	Ø 1.0 mm X 7 cm
DAG-103.13	Ø 1.5 mm X 7 cm
DAG-103.14	Ø 1.5 mm X 9 cm

Luer lock

DAG-104.11	Ø 1.0 mm X 4 cm
DAG-104.12	Ø 1.0 mm X 5 cm
DAG-104.13	Ø 1.0 mm X 7 cm
DAG-104.14	Ø 1.5 mm X 7 cm
DAG-104.15	Ø 1.5 mm X 9 cm
DAG-104.16	Ø 1.5 mm X 10 cm
DAG-104.17	Ø 1.5 mm X 12 cm
DAG-104.18	Ø 2.0 mm X 10 cm
DAG-104.19	Ø 2.0 mm X 12 cm
DAG-104.20	Ø 2.0 mm X 15 cm
DAG-104.21	Ø 2.0 mm X 20 cm
DAG-104.22	Ø 2.5 mm X 15 cm
DAG-104.23	Ø 2.5 mm X 20 cm
DAG-104.24	Ø 3.0 mm X 15 cm
DAG-104.25	Ø 3.0 mm X 20 cm


Luer lock (Stainless Steel)

DAG-105.21	Ø 1.0 mm X 4 cm
DAG-105.22	Ø 1.0 mm X 5 cm
DAG-105.23	Ø 1.0 mm X 7 cm
DAG-105.24	Ø 1.5 mm X 7 cm
DAG-105.25	Ø 1.5 mm X 9 cm
DAG-105.26	Ø 1.5 mm X 10 cm
DAG-105.27	Ø 1.5 mm X 12 cm
DAG-105.28	Ø 2.0 mm X 10 cm
DAG-105.29	Ø 2.0 mm X 12 cm
DAG-105.30	Ø 2.0 mm X 15 cm
DAG-105.31	Ø 2.0 mm X 20 cm
DAG-105.32	Ø 2.5 mm X 15 cm
DAG-105.33	Ø 2.5 mm X 20 cm
DAG-105.34	Ø 3.0 mm X 15 cm
DAG-105.35	Ø 3.0 mm X 20 cm

Power Cannula

DAG-106.31	Ø 2 mm X 15 cm
DAG-106.32	Ø 2 mm X 20 cm
DAG-106.33	Ø 2 mm X 25 cm
DAG-106.34	Ø 2 mm X 30 cm
DAG-106.35	Ø 2 mm X 35 cm
DAG-106.36	Ø 3 mm X 15 cm
DAG-106.37	Ø 3 mm X 20 cm
DAG-106.38	Ø 3 mm X 25 cm
DAG-106.39	Ø 3 mm X 30 cm
DAG-106.40	Ø 3 mm X 35 cm
DAG-106.41	Ø 4 mm X 15 cm
DAG-106.42	Ø 4 mm X 20 cm
DAG-106.43	Ø 4 mm X 25 cm
DAG-106.44	Ø 4 mm X 30 cm
DAG-106.45	Ø 4 mm X 35 cm
DAG-106.46	Ø 5 mm X 15 cm
DAG-106.47	Ø 5 mm X 20 cm
DAG-106.48	Ø 5 mm X 25 cm
DAG-106.49	Ø 5 mm X 30 cm
DAG-106.50	Ø 5 mm X 35 cm
DAG-106.51	Ø 6 mm X 15 cm
DAG-106.52	Ø 6 mm X 20 cm
DAG-106.53	Ø 6 mm X 25 cm
DAG-106.54	Ø 6 mm X 30 cm
DAG-106.55	Ø 6 mm X 35 cm

Infiltration Cannulas



Tumescent Infiltrator Angled

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60ccTommeyhub

DAG-101.21	Ø 2.0 mm X 15 cm
DAG-101.22	Ø 2.5 mm X 20 cm
DAG-101.23	Ø 3.0 mm X 20 cm
DAG-101.24	Ø 3.0 mm X 25 cm

Treaded Hub

DAG-102.21	Ø 2.0mm X 15 cm
DAG-102.22	Ø 3.0mm X 15 cm
DAG-102.23	Ø 3.0mm X 20 cm
DAG-102.24	Ø 3.0mm X 30 cm
DAG-102.25	Ø 4.0mm X 20 cm
DAG-102.26	Ø 4.0mm X 30 cm

Ultra Luer lock

DAG-103.21	Ø 1.0 mm X 5 cm
DAG-103.22	Ø 1.0 mm X 7 cm
DAG-103.23	Ø 1.5 mm X 7 cm
DAG-103.24	Ø 1.5 mm X 9 cm

Luer lock

DAG-104.31	Ø 1.0 mm X 4 cm
DAG-104.32	Ø 1.0 mm X 5 cm
DAG-104.33	Ø 1.0 mm X 7 cm
DAG-104.34	Ø 1.5 mm X 7 cm
DAG-104.35	Ø 1.5 mm X 9 cm
DAG-104.36	Ø 1.5 mm X 10 cm
DAG-104.37	Ø 1.5 mm X 12 cm
DAG-104.38	Ø 2.0 mm X 10 cm
DAG-104.39	Ø 2.0 mm X 12 cm
DAG-104.40	Ø 2.0 mm X 15 cm
DAG-104.41	Ø 2.0 mm X 20 cm
DAG-104.42	Ø 2.5 mm X 15 cm
DAG-104.43	Ø 2.5 mm X 20 cm
DAG-104.44	Ø 3.0 mm X 15 cm
DAG-104.45	Ø 3.0 mm X 20 cm


Luer lock (Stainless Steel)

DAG-105.41	Ø 1.0 mm X 4 cm
DAG-105.42	Ø 1.0 mm X 5 cm
DAG-105.43	Ø 1.0 mm X 7 cm
DAG-105.44	Ø 1.5 mm X 7 cm
DAG-105.45	Ø 1.5 mm X 9 cm
DAG-105.46	Ø 1.5 mm X 10 cm
DAG-105.47	Ø 1.5 mm X 12 cm
DAG-105.48	Ø 2.0 mm X 10 cm
DAG-105.49	Ø 2.0 mm X 12 cm
DAG-105.50	Ø 2.0 mm X 15 cm
DAG-105.51	Ø 2.0 mm X 20 cm
DAG-105.52	Ø 2.5 mm X 15 cm
DAG-105.53	Ø 2.5 mm X 20 cm
DAG-105.54	Ø 3.0 mm X 15 cm
DAG-105.55	Ø 3.0 mm X 20 cm

Power Cannula

DAG-106.61	Ø 2 mm X 15 cm
DAG-106.62	Ø 2 mm X 20 cm
DAG-106.63	Ø 2 mm X 25 cm
DAG-106.64	Ø 2 mm X 30 cm
DAG-106.65	Ø 2 mm X 35 cm
DAG-106.66	Ø 3 mm X 15 cm
DAG-106.67	Ø 3 mm X 20 cm
DAG-106.68	Ø 3 mm X 25 cm
DAG-106.69	Ø 3 mm X 30 cm
DAG-106.70	Ø 3 mm X 35 cm
DAG-106.71	Ø 4 mm X 15 cm
DAG-106.72	Ø 4 mm X 20 cm
DAG-106.73	Ø 4 mm X 25 cm
DAG-106.74	Ø 4 mm X 30 cm
DAG-106.75	Ø 4 mm X 35 cm
DAG-106.76	Ø 5 mm X 15 cm
DAG-106.77	Ø 5 mm X 20 cm
DAG-106.78	Ø 5 mm X 25 cm
DAG-106.79	Ø 5 mm X 30 cm
DAG-106.80	Ø 5 mm X 35 cm
DAG-106.81	Ø 6 mm X 15 cm
DAG-106.82	Ø 6 mm X 20 cm
DAG-106.83	Ø 6 mm X 25 cm
DAG-106.84	Ø 6 mm X 30 cm
DAG-106.85	Ø 6 mm X 35 cm

Infiltration Cannulas



Facelift Infiltration Cannula

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60ccTommeyhub

DAG-101.31	Ø 2.0 mm X 15 cm
DAG-101.32	Ø 2.5 mm X 20 cm
DAG-101.33	Ø 3.0 mm X 20 cm
DAG-101.34	Ø 3.0 mm X 25 cm

Treaded Hub

DAG-102.31	Ø 2.0mm X 15 cm
DAG-102.32	Ø 3.0mm X 15 cm
DAG-102.33	Ø 3.0mm X 20 cm
DAG-102.34	Ø 3.0mm X 30 cm
DAG-102.35	Ø 4.0mm X 20 cm
DAG-102.36	Ø 4.0mm X 30 cm

Ultra Luer lock

DAG-103.31	Ø 1.0 mm X 5 cm
DAG-103.32	Ø 1.0 mm X 7 cm
DAG-103.33	Ø 1.5 mm X 7 cm
DAG-103.34	Ø 1.5 mm X 9 cm

Luer lock

DAG-104.51	Ø 1.0 mm X 4 cm
DAG-104.52	Ø 1.0 mm X 5 cm
DAG-104.53	Ø 1.0 mm X 7 cm
DAG-104.54	Ø 1.5 mm X 7 cm
DAG-104.55	Ø 1.5 mm X 9 cm
DAG-104.56	Ø 1.5 mm X 10 cm
DAG-104.57	Ø 1.5 mm X 12 cm
DAG-104.58	Ø 2.0 mm X 10 cm
DAG-104.59	Ø 2.0 mm X 12 cm
DAG-104.50	Ø 2.0 mm X 15 cm
DAG-104.51	Ø 2.0 mm X 20 cm
DAG-104.52	Ø 2.5 mm X 15 cm
DAG-104.53	Ø 2.5 mm X 20 cm
DAG-104.54	Ø 3.0 mm X 15 cm
DAG-104.55	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

DAG-105.61	Ø 1.0 mm X 4 cm
DAG-105.62	Ø 1.0 mm X 5 cm
DAG-105.63	Ø 1.0 mm X 7 cm
DAG-105.64	Ø 1.5 mm X 7 cm
DAG-105.65	Ø 1.5 mm X 9 cm
DAG-105.66	Ø 1.5 mm X 10 cm
DAG-105.67	Ø 1.5 mm X 12 cm
DAG-105.68	Ø 2.0 mm X 10 cm
DAG-105.69	Ø 2.0 mm X 12 cm
DAG-105.70	Ø 2.0 mm X 15 cm
DAG-105.71	Ø 2.0 mm X 20 cm
DAG-105.72	Ø 2.5 mm X 15 cm
DAG-105.73	Ø 2.5 mm X 20 cm
DAG-105.74	Ø 3.0 mm X 15 cm
DAG-105.75	Ø 3.0 mm X 20 cm

Power Cannula

DAG-106.91	Ø 2 mm X 15 cm
DAG-106.92	Ø 2 mm X 20 cm
DAG-106.93	Ø 2 mm X 25 cm
DAG-106.94	Ø 2 mm X 30 cm
DAG-106.95	Ø 2 mm X 35 cm
DAG-106.96	Ø 3 mm X 15 cm
DAG-106.97	Ø 3 mm X 20 cm
DAG-106.98	Ø 3 mm X 25 cm
DAG-106.99	Ø 3 mm X 30 cm
DAG-106.100	Ø 3 mm X 35 cm
DAG-106.101	Ø 4 mm X 15 cm
DAG-106.102	Ø 4 mm X 20 cm
DAG-106.103	Ø 4 mm X 25 cm
DAG-106.104	Ø 4 mm X 30 cm
DAG-106.105	Ø 4 mm X 35 cm
DAG-106.106	Ø 5 mm X 15 cm
DAG-106.107	Ø 5 mm X 20 cm
DAG-106.108	Ø 5 mm X 25 cm
DAG-106.109	Ø 5 mm X 30 cm
DAG-106.110	Ø 5 mm X 35 cm
DAG-106.111	Ø 6 mm X 15 cm
DAG-106.112	Ø 6 mm X 20 cm
DAG-106.113	Ø 6 mm X 25 cm
DAG-106.114	Ø 6 mm X 30 cm
DAG-106.115	Ø 6 mm X 35 cm




SMART *life*
with



GENERAL SUCTION Cannulas

DAG SURGICAL CORPORATION General suction cannula are ideal tool for large volume Harvesting with ease perfection. These cannulas are designed to facilitate you to get perfectly designed goals in minimum time..



General Suction Cannulas



Stevens Style Cannula



60ccTommeyhub

DAG-201.01	Ø 2.0 mm X 15 cm
DAG-201.02	Ø 2.5 mm X 20 cm
DAG-201.03	Ø 3.0 mm X 20 cm
DAG-201.04	Ø 3.0 mm X 25 cm

Treaded Hub

DAG-202.01	Ø 3.0 mm X 20 cm
DAG-202.02	Ø 3.0 mm X 30 cm
DAG-202.03	Ø 4.0 mm X 30 cm
DAG-202.04	Ø 5.0 mm X 30 cm
DAG-202.05	Ø 6.0 mm X 30 cm

Ultra Luer lock

DAG-203.01	Ø 1.0 mm X 5 cm
DAG-203.02	Ø 1.0 mm X 7 cm
DAG-203.03	Ø 1.5 mm X 7 cm
DAG-203.04	Ø 1.5 mm X 9 cm
DAG-203.05	Ø 1.5 mm X 7 cm

Luer lock

DAG-204.01	Ø 1.0 mm X 4 cm
DAG-204.02	Ø 1.0 mm X 5 cm
DAG-204.03	Ø 1.0 mm X 7 cm
DAG-204.04	Ø 1.5 mm X 7 cm
DAG-204.05	Ø 1.5 mm X 9 cm
DAG-204.06	Ø 1.5 mm X 10 cm
DAG-204.07	Ø 1.5 mm X 12 cm
DAG-204.08	Ø 2.0 mm X 10 cm
DAG-204.09	Ø 2.0 mm X 12 cm
DAG-204.10	Ø 2.0 mm X 15 cm
DAG-204.11	Ø 2.0 mm X 20 cm
DAG-204.12	Ø 2.5 mm X 15 cm
DAG-204.13	Ø 2.5 mm X 20 cm
DAG-204.14	Ø 3.0 mm X 15 cm
DAG-204.15	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

DAG-205.01	Ø 1.0 mm X 4 cm
DAG-205.02	Ø 1.0 mm X 5 cm
DAG-205.03	Ø 1.0 mm X 7 cm
DAG-205.04	Ø 1.5 mm X 7 cm
DAG-205.05	Ø 1.5 mm X 9 cm
DAG-205.06	Ø 1.5 mm X 10 cm
DAG-205.07	Ø 1.5 mm X 12 cm
DAG-205.08	Ø 2.0 mm X 10 cm
DAG-205.09	Ø 2.0 mm X 12 cm
DAG-205.10	Ø 2.0 mm X 15 cm
DAG-205.11	Ø 2.0 mm X 20 cm
DAG-205.12	Ø 2.5 mm X 15 cm
DAG-205.13	Ø 2.5 mm X 20 cm
DAG-205.14	Ø 3.0 mm X 15 cm
DAG-205.15	Ø 3.0 mm X 20 cm

Power Cannula


DAG-206.01	Ø 2 mm X 15 cm
DAG-206.02	Ø 2 mm X 20 cm
DAG-206.03	Ø 2 mm X 25 cm
DAG-206.04	Ø 2 mm X 30 cm
DAG-206.05	Ø 2 mm X 35 cm
DAG-206.06	Ø 3 mm X 15 cm
DAG-206.07	Ø 3 mm X 20 cm
DAG-206.08	Ø 3 mm X 25 cm
DAG-206.09	Ø 3 mm X 30 cm
DAG-206.10	Ø 3 mm X 35 cm
DAG-206.11	Ø 4 mm X 15 cm
DAG-206.12	Ø 4 mm X 20 cm
DAG-206.13	Ø 4 mm X 25 cm
DAG-206.14	Ø 4 mm X 30 cm
DAG-206.15	Ø 4 mm X 35 cm
DAG-206.16	Ø 5 mm X 15 cm
DAG-206.17	Ø 5 mm X 20 cm
DAG-206.18	Ø 5 mm X 25 cm
DAG-206.19	Ø 5 mm X 30 cm
DAG-206.20	Ø 5 mm X 35 cm
DAG-206.21	Ø 6 mm X 15 cm
DAG-206.22	Ø 6 mm X 20 cm
DAG-206.23	Ø 6 mm X 25 cm
DAG-206.24	Ø 6 mm X 30 cm
DAG-206.25	Ø 6 mm X 35 cm

General Suction Cannulas



Cobra Bibevel Style

More aggressive than the Standard, Pyramid or Mercedes, this cannula has a tapered tip with two ports and works effectively within dense fibrous tissue.



60ccTommeyhub

DAG-201.11	Ø 2.5 mm X 15 cm5
DAG-201.12	Ø 2.5 mm X 20 cm
DAG-201.13	Ø 3.0 mm X 15 cm
DAG-201.14	Ø 3.0 mm X 20 cm
DAG-201.15	Ø 3.0 mm X 25 cm
DAG-201.16	Ø 4.0 mm X 25 cm
DAG-201.17	Ø 4.0 mm X 30 cm
DAG-201.18	Ø 4.0 mm X 35 cm
DAG-201.19	Ø 4.5 mm X 30 cm
DAG-201.20	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.11	Ø 3.0 mm X 15 cm
DAG-202.12	Ø 3.0 mm X 25 cm
DAG-202.13	Ø 4.0 mm X 15 cm
DAG-202.14	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.11	Ø 2.0 mm X 15 cm
DAG-203.12	Ø 2.5 mm X 15 cm
DAG-203.13	Ø 2.5 mm X 20 cm
DAG-203.14	Ø 3.0 mm X 15 cm
DAG-203.15	Ø 3.0 mm X 20 cm
DAG-203.16	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.21	Ø 2.0 mm X 15 cm
DAG-204.22	Ø 2.5 mm X 15 cm
DAG-204.23	Ø 2.5 mm X 20 cm
DAG-204.24	Ø 3.0 mm X 15 cm
DAG-204.25	Ø 3.0 mm X 20 cm
DAG-204.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

DAG-205.21	Ø 2.0 mm X 15 cm
DAG-205.22	Ø 2.5 mm X 15 cm
DAG-205.23	Ø 2.5 mm X 20 cm
DAG-205.24	Ø 3.0 mm X 15 cm
DAG-205.25	Ø 3.0 mm X 20 cm
DAG-205.26	Ø 3.0 mm X 25 cm

Power Cannula


DAG-206.31	Ø 2mm X 15 cm
DAG-206.32	Ø 2mm X 20 cm
DAG-206.33	Ø 2mm X 25 cm
DAG-206.34	Ø 2mm X 30 cm
DAG-206.35	Ø 2mm X 35 cm
DAG-206.36	Ø 3mm X 15 cm
DAG-206.37	Ø 3mm X 20 cm
DAG-206.38	Ø 3mm X 25 cm
DAG-206.39	Ø 3mm X 30 cm
DAG-206.40	Ø 3mm X 35 cm
DAG-206.41	Ø 4mm X 15 cm
DAG-206.42	Ø 4mm X 20 cm
DAG-206.43	Ø 4mm X 25 cm
DAG-206.44	Ø 4mm X 30 cm
DAG-206.45	Ø 4mm X 35 cm
DAG-206.46	Ø 5mm X 15 cm
DAG-206.47	Ø 5mm X 20 cm
DAG-206.48	Ø 5mm X 25 cm
DAG-206.49	Ø 5mm X 30 cm
DAG-206.50	Ø 5mm X 35 cm
DAG-206.51	Ø 6mm X 15 cm
DAG-206.52	Ø 6mm X 20 cm
DAG-206.53	Ø 6mm X 25 cm
DAG-206.54	Ø 6mm X 30 cm
DAG-206.55	Ø 6mm X 35 cm

General Suction Cannulas



Cobra Style Round Tip

Less aggressive than the Cobra Bibevel this cell friendly cannula has a rounded tip with two ports and works effectively within dense fibrous tissue during liposuction procedures.



60ccTommeyhub

DAG-201.21	Ø 2.5 mm X 15 cm
DAG-201.22	Ø 2.5 mm X 20 cm
DAG-201.23	Ø 3.0 mm X 15 cm
DAG-201.24	Ø 3.0 mm X 20 cm
DAG-201.25	Ø 3.0 mm X 25 cm
DAG-201.26	Ø 4.0 mm X 25 cm
DAG-201.27	Ø 4.0 mm X 30 cm
DAG-201.28	Ø 4.0 mm X 35 cm
DAG-201.29	Ø 4.5 mm X 30 cm
DAG-201.30	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.21	Ø 3.0 mm X 15 cm
DAG-202.22	Ø 3.0 mm X 25 cm
DAG-202.23	Ø 4.0 mm X 15 cm
DAG-202.24	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.21	Ø 2.0 mm X 15 cm
DAG-203.22	Ø 2.5 mm X 15 cm
DAG-203.23	Ø 2.5 mm X 20 cm
DAG-203.24	Ø 3.0 mm X 15 cm
DAG-203.25	Ø 3.0 mm X 20 cm
DAG-203.26	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.31	Ø 2.0 mm X 15 cm
DAG-204.32	Ø 2.5 mm X 15 cm
DAG-204.33	Ø 2.5 mm X 20 cm
DAG-204.34	Ø 3.0 mm X 15 cm
DAG-204.35	Ø 3.0 mm X 20 cm
DAG-204.36	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

DAG-205.31	Ø 2.0 mm X 15 cm
DAG-205.32	Ø 2.5 mm X 15 cm
DAG-205.33	Ø 2.5 mm X 20 cm
DAG-205.34	Ø 3.0 mm X 15 cm
DAG-205.35	Ø 3.0 mm X 20 cm
DAG-205.36	Ø 3.0 mm X 25 cm

Power Cannula


DAG-206.61	Ø 2 mm X 15 cm
DAG-206.62	Ø 2 mm X 20 cm
DAG-206.63	Ø 2 mm X 25 cm
DAG-206.64	Ø 2 mm X 30 cm
DAG-206.65	Ø 2 mm X 35 cm
DAG-206.66	Ø 3 mm X 15 cm
DAG-206.67	Ø 3 mm X 20 cm
DAG-206.68	Ø 3 mm X 25 cm
DAG-206.69	Ø 3 mm X 30 cm
DAG-206.70	Ø 3 mm X 35 cm
DAG-206.71	Ø 4 mm X 15 cm
DAG-206.72	Ø 4 mm X 20 cm
DAG-206.73	Ø 4 mm X 25 cm
DAG-206.74	Ø 4 mm X 30 cm
DAG-206.75	Ø 4 mm X 35 cm
DAG-206.76	Ø 5 mm X 15 cm
DAG-206.77	Ø 5 mm X 20 cm
DAG-206.78	Ø 5 mm X 25 cm
DAG-206.79	Ø 5 mm X 30 cm
DAG-206.80	Ø 5 mm X 35 cm
DAG-206.81	Ø 6 mm X 15 cm
DAG-206.82	Ø 6 mm X 20 cm
DAG-206.83	Ø 6 mm X 25 cm
DAG-206.84	Ø 6 mm X 30 cm
DAG-206.85	Ø 6 mm X 35 cm

General Suction Cannulas



Standard, Single Hole

This is the most common cell friendly cannula design for soft tissue removal. The single port is located 180 degrees from the thumb rest on the hub.



60ccTommeyhub

DAG-201.31	Ø 2.5 mm X 15 cm
DAG-201.32	Ø 2.5 mm X 20 cm
DAG-201.33	Ø 3.0 mm X 15 cm
DAG-201.34	Ø 3.0 mm X 20 cm
DAG-201.35	Ø 3.0 mm X 25 cm
DAG-201.36	Ø 4.0 mm X 25 cm
DAG-201.37	Ø 4.0 mm X 30 cm
DAG-201.38	Ø 4.0 mm X 35 cm
DAG-201.39	Ø 4.5 mm X 30 cm
DAG-201.40	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.31	Ø 3.0 mm X 15 cm
DAG-202.32	Ø 3.0 mm X 25 cm
DAG-202.33	Ø 4.0 mm X 15 cm
DAG-202.34	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.31	Ø 1.5 mm X 12 cm
DAG-203.32	Ø 2.0 mm X 15 cm
DAG-203.33	Ø 2.0 mm X 20 cm
DAG-203.34	Ø 2.5 mm X 15 cm
DAG-203.35	Ø 2.5 mm X 20 cm
DAG-203.36	Ø 3.0 mm X 15 cm
DAG-203.37	Ø 3.0 mm X 20 cm
DAG-203.38	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.41	Ø 1.5 mm X 12 cm
DAG-204.42	Ø 2.0 mm X 15 cm
DAG-204.43	Ø 2.0 mm X 20 cm
DAG-204.44	Ø 2.5 mm X 15 cm
DAG-204.45	Ø 2.5 mm X 20 cm
DAG-204.46	Ø 3.0 mm X 15 cm
DAG-204.47	Ø 3.0 mm X 20 cm
DAG-204.48	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-205.41	Ø 1.5 mm X 12 cm
DAG-205.42	Ø 2.0 mm X 15 cm
DAG-205.43	Ø 2.0 mm X 20 cm
DAG-205.44	Ø 2.5 mm X 15 cm
DAG-205.45	Ø 2.5 mm X 20 cm
DAG-205.46	Ø 3.0 mm X 15 cm
DAG-205.47	Ø 3.0 mm X 20 cm
DAG-205.48	Ø 3.0 mm X 25 cm

Power Cannula

DAG-206.91	Ø 2 mm X 15 cm
DAG-206.92	Ø 2 mm X 20 cm
DAG-206.93	Ø 2 mm X 25 cm
DAG-206.94	Ø 2 mm X 30 cm
DAG-206.95	Ø 2 mm X 35 cm
DAG-206.96	Ø 3 mm X 15 cm
DAG-206.97	Ø 3 mm X 20 cm
DAG-206.98	Ø 3 mm X 25 cm
DAG-206.99	Ø 3 mm X 30 cm
DAG-206.100	Ø 3 mm X 35 cm
DAG-206.101	Ø 4 mm X 15 cm
DAG-206.102	Ø 4 mm X 20 cm
DAG-206.103	Ø 4 mm X 25 cm
DAG-206.104	Ø 4 mm X 30 cm
DAG-206.105	Ø 4 mm X 35 cm
DAG-206.106	Ø 5 mm X 15 cm
DAG-206.107	Ø 5 mm X 20 cm
DAG-206.108	Ø 5 mm X 25 cm
DAG-206.109	Ø 5 mm X 30 cm
DAG-206.110	Ø 5 mm X 35 cm
DAG-206.111	Ø 6 mm X 15 cm
DAG-206.112	Ø 6 mm X 20 cm
DAG-206.113	Ø 6 mm X 25 cm
DAG-206.114	Ø 6 mm X 30 cm
DAG-206.115	Ø 6 mm X 35 cm

General Suction Cannulas



Mercedes

This circumferential three port configuration cell friendly liposuction cannula is effective for multi-plane soft tissue harvesting.

60ccTommeyhub

DAG-201.51	Ø 2.5 mm X 15 cm
DAG-201.52	Ø 2.5 mm X 20 cm
DAG-201.53	Ø 3.0 mm X 15 cm
DAG-201.54	Ø 3.0 mm X 20 cm
DAG-201.55	Ø 3.0 mm X 25 cm
DAG-201.56	Ø 4.0 mm X 25 cm
DAG-201.57	Ø 4.0 mm X 30 cm
DAG-201.58	Ø 4.0 mm X 35 cm
DAG-201.59	Ø 4.5 mm X 30 cm
DAG-201.60	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.51	Ø 3.0 mm X 15 cm
DAG-202.52	Ø 3.0 mm X 25 cm
DAG-202.53	Ø 4.0 mm X 15 cm
DAG-202.54	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.51	Ø 1.5 mm X 12 cm
DAG-203.52	Ø 2.0 mm X 15 cm
DAG-203.53	Ø 2.5 mm X 15 cm
DAG-203.54	Ø 2.5 mm X 20 cm
DAG-203.55	Ø 3.0 mm X 15 cm
DAG-203.56	Ø 3.0 mm X 20 cm
DAG-203.57	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.61	Ø 1.5 mm X 12 cm
DAG-204.62	Ø 2.0 mm X 15 cm
DAG-204.63	Ø 2.5 mm X 15 cm
DAG-204.64	Ø 2.5 mm X 20 cm
DAG-204.65	Ø 3.0 mm X 15 cm
DAG-204.66	Ø 3.0 mm X 20 cm
DAG-204.67	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-205.61	Ø 1.5 mm X 12 cm
DAG-205.62	Ø 2.0 mm X 15 cm
DAG-205.63	Ø 2.0 mm X 20 cm
DAG-205.64	Ø 2.5 mm X 15 cm
DAG-205.65	Ø 2.5 mm X 20 cm
DAG-205.66	Ø 3.0 mm X 15 cm
DAG-205.67	Ø 3.0 mm X 20 cm
DAG-205.68	Ø 3.0 mm X 25 cm

Power Cannula

DAG-206.151	Ø 2 mm X 15 cm
DAG-206.152	Ø 2 mm X 20 cm
DAG-206.153	Ø 2 mm X 25 cm
DAG-206.154	Ø 2 mm X 30 cm
DAG-206.155	Ø 2 mm X 35 cm
DAG-206.156	Ø 3 mm X 15 cm
DAG-206.157	Ø 3 mm X 20 cm
DAG-206.158	Ø 3 mm X 25 cm
DAG-206.159	Ø 3 mm X 30 cm
DAG-206.160	Ø 3 mm X 35 cm
DAG-206.161	Ø 4 mm X 15 cm
DAG-206.162	Ø 4 mm X 20 cm
DAG-206.163	Ø 4 mm X 25 cm
DAG-206.164	Ø 4 mm X 30 cm
DAG-206.165	Ø 4 mm X 35 cm
DAG-206.166	Ø 5 mm X 15 cm
DAG-206.167	Ø 5 mm X 20 cm
DAG-206.168	Ø 5 mm X 25 cm
DAG-206.169	Ø 5 mm X 30 cm
DAG-206.170	Ø 5 mm X 35 cm
DAG-206.171	Ø 6 mm X 15 cm
DAG-206.172	Ø 6 mm X 20 cm
DAG-206.173	Ø 6 mm X 25 cm
DAG-206.174	Ø 6 mm X 30 cm
DAG-206.175	Ø 6 mm X 35 cm

General Suction Cannulas



Pyramid Holes Cannula

This commonly used design has three ports located 180 degrees from the thumb print. It is a slightly more aggressive cannula than the Standard during liposuction procedures.

60ccTommeyhub

DAG-201.41	Ø 2.5 mm X 15 cm
DAG-201.42	Ø 2.5 mm X 20 cm
DAG-201.43	Ø 3.0 mm X 15 cm
DAG-201.44	Ø 3.0 mm X 20 cm
DAG-201.45	Ø 3.0 mm X 25 cm
DAG-201.46	Ø 4.0 mm X 25 cm
DAG-201.47	Ø 4.0 mm X 30 cm
DAG-201.48	Ø 4.0 mm X 35 cm
DAG-201.49	Ø 4.5 mm X 30 cm
DAG-201.50	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.41	Ø 3.0 mm X 15 cm
DAG-202.42	Ø 3.0 mm X 25 cm
DAG-202.43	Ø 4.0 mm X 15 cm
DAG-202.44	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.41	Ø 1.5 mm X 12 cm
DAG-203.42	Ø 2.0 mm X 15 cm
DAG-203.43	Ø 2.5 mm X 15 cm
DAG-203.44	Ø 2.5 mm X 20 cm
DAG-203.45	Ø 3.0 mm X 15 cm
DAG-203.46	Ø 3.0 mm X 20 cm
DAG-203.47	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.51	Ø 1.5 mm X 12 cm
DAG-204.52	Ø 2.0 mm X 15 cm
DAG-204.53	Ø 2.5 mm X 15 cm
DAG-204.54	Ø 2.5 mm X 20 cm
DAG-204.55	Ø 3.0 mm X 15 cm
DAG-204.56	Ø 3.0 mm X 20 cm
DAG-204.57	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

DAG-205.51	Ø 1.5 mm X 12 cm
DAG-205.52	Ø 2.0 mm X 15 cm
DAG-205.53	Ø 2.0 mm X 20 cm
DAG-205.54	Ø 2.5 mm X 15 cm
DAG-205.55	Ø 2.5 mm X 20 cm
DAG-205.56	Ø 3.0 mm X 15 cm
DAG-205.57	Ø 3.0 mm X 20 cm
DAG-205.58	Ø 3.0 mm X 25 cm

Power Cannula

DAG-206.121	Ø 2 mm X 15 cm
DAG-206.122	Ø 2 mm X 20 cm
DAG-206.123	Ø 2 mm X 25 cm
DAG-206.124	Ø 2 mm X 30 cm
DAG-206.125	Ø 2 mm X 35 cm
DAG-206.126	Ø 3 mm X 15 cm
DAG-206.127	Ø 3 mm X 20 cm
DAG-206.128	Ø 3 mm X 25 cm
DAG-206.129	Ø 3 mm X 30 cm
DAG-206.130	Ø 3 mm X 35 cm
DAG-206.131	Ø 4 mm X 15 cm
DAG-206.132	Ø 4 mm X 20 cm
DAG-206.133	Ø 4 mm X 25 cm
DAG-206.134	Ø 4 mm X 30 cm
DAG-206.135	Ø 4 mm X 35 cm
DAG-206.136	Ø 5 mm X 15 cm
DAG-206.137	Ø 5 mm X 20 cm
DAG-206.138	Ø 5 mm X 25 cm
DAG-206.139	Ø 5 mm X 30 cm
DAG-206.140	Ø 5 mm X 35 cm
DAG-206.141	Ø 6 mm X 15 cm
DAG-206.142	Ø 6 mm X 20 cm
DAG-206.143	Ø 6 mm X 25 cm
DAG-206.144	Ø 6 mm X 30 cm
DAG-206.145	Ø 6 mm X 35 cm

General Suction Cannulas



3 Spiral Port Cannula (Carraway Style)

Spiral three port design on this Cell Friendly liposuction cannula is effective for multi-plane adipose tissue harvesting.

60ccTommyhub

DAG-201.61	Ø 2.5 mm X 15 cm
DAG-201.62	Ø 2.5 mm X 20 cm
DAG-201.63	Ø 3.0 mm X 15 cm
DAG-201.64	Ø 3.0 mm X 20 cm
DAG-201.65	Ø 3.0 mm X 25 cm
DAG-201.66	Ø 4.0 mm X 25 cm
DAG-201.67	Ø 4.0 mm X 30 cm
DAG-201.68	Ø 4.0 mm X 35 cm
DAG-201.69	Ø 4.5 mm X 30 cm
DAG-201.70	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.61	Ø 3.0 mm X 15 cm
DAG-202.62	Ø 3.0 mm X 25 cm
DAG-202.63	Ø 4.0 mm X 15 cm
DAG-202.64	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.61	Ø 1.5 mm X 12 cm
DAG-203.62	Ø 2.0 mm X 15 cm
DAG-203.63	Ø 2.5 mm X 15 cm
DAG-203.64	Ø 2.5 mm X 20 cm
DAG-203.65	Ø 3.0 mm X 15 cm
DAG-203.66	Ø 3.0 mm X 20 cm
DAG-203.67	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.71	Ø 1.5 mm X 12 cm
DAG-204.72	Ø 2.0 mm X 15 cm
DAG-204.73	Ø 2.5 mm X 15 cm
DAG-204.74	Ø 2.5 mm X 20 cm
DAG-204.75	Ø 3.0 mm X 15 cm
DAG-204.76	Ø 3.0 mm X 20 cm
DAG-204.77	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-205.71	Ø 1.5 mm X 12 cm
DAG-205.72	Ø 2.0 mm X 15 cm
DAG-205.73	Ø 2.0 mm X 20 cm
DAG-205.74	Ø 2.5 mm X 15 cm
DAG-205.75	Ø 2.5 mm X 20 cm
DAG-205.76	Ø 3.0 mm X 15 cm
DAG-205.77	Ø 3.0 mm X 20 cm
DAG-205.78	Ø 3.0 mm X 25 cm

Power Cannula

DAG-206.181	Ø 2 mm X 15 cm
DAG-206.182	Ø 2 mm X 20 cm
DAG-206.183	Ø 2 mm X 25 cm
DAG-206.184	Ø 2 mm X 30 cm
DAG-206.185	Ø 2 mm X 35 cm
DAG-206.186	Ø 3 mm X 15 cm
DAG-206.187	Ø 3 mm X 20 cm
DAG-206.188	Ø 3 mm X 25 cm
DAG-206.189	Ø 3 mm X 30 cm
DAG-206.190	Ø 3 mm X 35 cm
DAG-206.191	Ø 4 mm X 15 cm
DAG-206.192	Ø 4 mm X 20 cm
DAG-206.193	Ø 4 mm X 25 cm
DAG-206.194	Ø 4 mm X 30 cm
DAG-206.195	Ø 4 mm X 35 cm
DAG-206.196	Ø 5 mm X 15 cm
DAG-206.197	Ø 5 mm X 20 cm
DAG-206.198	Ø 5 mm X 25 cm
DAG-206.199	Ø 5 mm X 30 cm
DAG-206.200	Ø 5 mm X 35 cm
DAG-206.201	Ø 6 mm X 15 cm
DAG-206.202	Ø 6 mm X 20 cm
DAG-206.203	Ø 6 mm X 25 cm
DAG-206.204	Ø 6 mm X 30 cm
DAG-206.205	Ø 6 mm X 35 cm

General Suction Cannulas



Six Spiral Port (Adjust below 3 Spiral Ports)

Six port spiral configuration on this Cell Friendly liposuction cannula harvests maximum tissue volume in minimal time.

60ccsmarthub

DAG-201.91	Ø 2.5 mm X 15 cm
DAG-201.92	Ø 2.5 mm X 20 cm
DAG-201.93	Ø 3.0 mm X 15 cm
DAG-201.94	Ø 3.0 mm X 25 cm
DAG-201.95	Ø 4.0 mm X 25 cm

Treaded Hub

DAG-202.81	Ø 3.0 mm X 15 cm
DAG-202.82	Ø 3.0 mm X 25 cm
DAG-202.83	Ø 4.0 mm X 15 cm
DAG-202.84	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.81	Ø 2.5 mm X 15 cm
DAG-203.82	Ø 3.0 mm X 15 cm
DAG-203.83	Ø 3.0 mm X 25 cm
DAG-203.84	Ø 4.0 mm X 15 cm
DAG-203.85	Ø 4.0 mm X 25 cm

Luer lock

DAG-204.91	Ø 2.5 mm X 15 cm
DAG-204.92	Ø 3.0 mm X 15 cm
DAG-204.93	Ø 3.0 mm X 25 cm
DAG-204.94	Ø 4.0 mm X 15 cm
DAG-204.95	Ø 4.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-205.231	Ø 2.5 mm X 15 cm
DAG-205.232	Ø 3.0 mm X 15 cm
DAG-205.233	Ø 3.0 mm X 25 cm
DAG-205.234	Ø 4.0 mm X 15 cm
DAG-205.235	Ø 4.0 mm X 25 cm

Power Cannula

DAG-206.901	2mm OD X 15 cm
DAG-206.902	2mm OD X 20 cm
DAG-206.903	2mm OD X 25 cm
DAG-206.904	2mm OD X 30 cm
DAG-206.905	2mm OD X 35 cm
DAG-206.906	3mm OD X 15 cm
DAG-206.907	3mm OD X 20 cm
DAG-206.908	3mm OD X 25 cm
DAG-206.909	3mm OD X 30 cm
DAG-206.910	3mm OD X 35 cm
DAG-206.911	4mm OD X 15 cm
DAG-206.912	4mm OD X 20 cm
DAG-206.913	4mm OD X 25 cm
DAG-206.914	4mm OD X 30 cm
DAG-206.915	4mm OD X 35 cm
DAG-206.916	5mm OD X 15 cm
DAG-206.917	5mm OD X 20 cm
DAG-206.918	5mm OD X 25 cm
DAG-206.919	5mm OD X 30 cm
DAG-206.920	5mm OD X 35 cm
DAG-206.921	6mm OD X 15 cm
DAG-206.922	6mm OD X 20 cm
DAG-206.923	6mm OD X 25 cm
DAG-206.924	6mm OD X 30 cm
DAG-206.925	6mm OD X 35 cm

General Suction Cannulas



Three Port Cannula

This body sculpting instrument is excellent for harvesting small aliquots of tissue during liposuction procedures.

60ccTommeyhub

DAG-201.71	Ø 2.5 mm X 15 cm
DAG-201.72	Ø 2.5 mm X 20 cm
DAG-201.73	Ø 3.0 mm X 15 cm
DAG-201.74	Ø 3.0 mm X 20 cm
DAG-201.75	Ø 3.0 mm X 25 cm
DAG-201.76	Ø 4.0 mm X 25 cm
DAG-201.77	Ø 4.0 mm X 30 cm
DAG-201.78	Ø 4.0 mm X 35 cm
DAG-201.79	Ø 4.5 mm X 30 cm
DAG-201.80	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-202.71	Ø 3.0 mm X 15 cm
DAG-202.72	Ø 3.0 mm X 25 cm
DAG-202.73	Ø 4.0 mm X 15 cm
DAG-202.74	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-203.71	Ø 1.5 mm X 12 cm
DAG-203.72	Ø 2.0 mm X 15 cm
DAG-203.73	Ø 2.5 mm X 15 cm
DAG-203.74	Ø 2.5 mm X 20 cm
DAG-203.75	Ø 3.0 mm X 15 cm
DAG-203.76	Ø 3.0 mm X 20 cm
DAG-203.77	Ø 3.0 mm X 25 cm

Luer lock

DAG-204.81	Ø 1.5 mm X 12 cm
DAG-204.82	Ø 2.0 mm X 15 cm
DAG-204.83	Ø 2.5 mm X 15 cm
DAG-204.84	Ø 2.5 mm X 20 cm
DAG-204.85	Ø 3.0 mm X 15 cm
DAG-204.86	Ø 3.0 mm X 20 cm
DAG-204.87	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-205.81	Ø 1.5 mm X 12 cm
DAG-205.82	Ø 2.0 mm X 15 cm
DAG-205.83	Ø 2.0 mm X 20 cm
DAG-205.84	Ø 2.5 mm X 15 cm
DAG-205.85	Ø 2.5 mm X 20 cm
DAG-205.86	Ø 3.0 mm X 15 cm
DAG-205.87	Ø 3.0 mm X 20 cm
DAG-205.88	Ø 3.0 mm X 25 cm

Power Cannula

DAG-206.211	Ø 2 mm X 15 cm
DAG-206.212	Ø 2 mm X 20 cm
DAG-206.213	Ø 2 mm X 25 cm
DAG-206.214	Ø 2 mm X 30 cm
DAG-206.215	Ø 2 mm X 35 cm
DAG-206.216	Ø 3 mm X 15 cm
DAG-206.217	Ø 3 mm X 20 cm
DAG-206.218	Ø 3 mm X 25 cm
DAG-206.219	Ø 3 mm X 30 cm
DAG-206.220	Ø 3 mm X 35 cm
DAG-206.221	Ø 4 mm X 15 cm
DAG-206.222	Ø 4 mm X 20 cm
DAG-206.223	Ø 4 mm X 25 cm
DAG-206.224	Ø 4 mm X 30 cm
DAG-206.225	Ø 4 mm X 35 cm
DAG-206.226	Ø 5 mm X 15 cm
DAG-206.227	Ø 5 mm X 20 cm
DAG-206.228	Ø 5 mm X 25 cm
DAG-206.229	Ø 5 mm X 30 cm
DAG-206.230	Ø 5 mm X 35 cm
DAG-206.231	Ø 6 mm X 15 cm
DAG-206.232	Ø 6 mm X 20 cm
DAG-206.233	Ø 6 mm X 25 cm
DAG-206.234	Ø 6 mm X 30 cm
DAG-206.235	Ø 6 mm X 35 cm

General Suction Cannulas



Basket Shaped Cannula

The Cell Friendly Basket Cannula is a versatile liposculpting tool. It can operate as an equalizer without vacuum to treat adipose tissue and skin irregularities in revision liposuction and fat transfer procedures. The Basket Cannula is also adept under vacuum as a fat harvesting cannula. For finishing, this cannula can be used as an injector for equalizing and adding adipose tissue to sculpt smoother contours.

60ccTommeyhub

DAG-201.81	Ø 3.0 mm X 15 cm
DAG-201.82	Ø 3.0 mm X 20 cm
DAG-201.83	Ø 3.0 mm X 25 cm
DAG-201.84	Ø 4.0 mm X 25 cm
DAG-201.85	Ø 4.0 mm X 30 cm
DAG-201.86	Ø 4.0 mm X 35 cm
DAG-201.87	Ø 5.0 mm X 25 cm
DAG-201.88	Ø 5.0 mm X 30 cm
DAG-201.89	Ø 5.0 mm X 35 cm
DAG-201.90	Ø 6.0 mm X 25 cm
DAG-201.91	Ø 6.0 mm X 30 cm
DAG-201.92	Ø 6.0 mm X 35 cm

Hole on Tip-Style I

Power Cannula

DAG-206.251	Ø 2 mm X 15 cm
DAG-206.252	Ø 2 mm X 20 cm
DAG-206.253	Ø 2 mm X 25 cm
DAG-206.254	Ø 2 mm X 30 cm
DAG-206.255	Ø 2 mm X 35 cm
DAG-206.256	Ø 3 mm X 15 cm
DAG-206.257	Ø 3 mm X 20 cm
DAG-206.258	Ø 3 mm X 25 cm
DAG-206.259	Ø 3 mm X 30 cm

DAG-206.260	Ø 3 mm X 35 cm
DAG-206.261	Ø 4 mm X 15 cm
DAG-206.262	Ø 4 mm X 20 cm
DAG-206.263	Ø 4 mm X 25 cm
DAG-206.264	Ø 4 mm X 30 cm
DAG-206.265	Ø 4 mm X 35 cm
DAG-206.266	Ø 5 mm X 15 cm
DAG-206.267	Ø 5 mm X 20 cm

DAG-206.268	Ø 5 mm X 25 cm
DAG-206.269	Ø 5 mm X 30 cm
DAG-206.270	Ø 5 mm X 35 cm
DAG-206.271	Ø 6 mm X 15 cm
DAG-206.272	Ø 6 mm X 20 cm
DAG-206.273	Ø 6 mm X 25 cm
DAG-206.274	Ø 6 mm X 30 cm
DAG-206.275	Ø 6 mm X 35 cm

Hole on Tip-Style II

Power Cannula

DAG-206.281	Ø 2 mm X 15 cm
DAG-206.282	Ø 2 mm X 20 cm
DAG-206.283	Ø 2 mm X 25 cm
DAG-206.284	Ø 2 mm X 30 cm
DAG-206.285	Ø 2 mm X 35 cm
DAG-206.286	Ø 3 mm X 15 cm
DAG-206.287	Ø 3 mm X 20 cm
DAG-206.288	Ø 3 mm X 25 cm
DAG-206.289	Ø 3 mm X 30 cm

DAG-206.290	Ø 3 mm X 35 cm
DAG-206.291	Ø 4 mm X 15 cm
DAG-206.292	Ø 4 mm X 20 cm
DAG-206.293	Ø 4 mm X 25 cm
DAG-206.294	Ø 4 mm X 30 cm
DAG-206.295	Ø 4 mm X 35 cm
DAG-206.296	Ø 5 mm X 15 cm
DAG-206.297	Ø 5 mm X 20 cm

DAG-206.298	Ø 5 mm X 25 cm
DAG-206.299	Ø 5 mm X 30 cm
DAG-206.300	Ø 5 mm X 35 cm
DAG-206.301	Ø 6 mm X 15 cm
DAG-206.302	Ø 6 mm X 20 cm
DAG-206.303	Ø 6 mm X 25 cm
DAG-206.304	Ø 6 mm X 30 cm
DAG-206.305	Ø 6 mm X 35 cm

General Suction Cannulas



Deep Sharp Cut Dissector

Power Cannula

DAG-206.311	Ø 2 mm X 15 cm	DAG-206.320	Ø 3 mm X 35 cm	DAG-206.328	Ø 5 mm X 25 cm
DAG-206.312	Ø 2 mm X 20 cm	DAG-206.321	Ø 4 mm X 15 cm	DAG-206.329	Ø 5 mm X 30 cm
DAG-206.313	Ø 2 mm X 25 cm	DAG-206.322	Ø 4 mm X 20 cm	DAG-206.330	Ø 5 mm X 35 cm
DAG-206.314	Ø 2 mm X 30 cm	DAG-206.323	Ø 4 mm X 25 cm	DAG-206.331	Ø 6 mm X 15 cm
DAG-206.315	Ø 2 mm X 35 cm	DAG-206.324	Ø 4 mm X 30 cm	DAG-206.332	Ø 6 mm X 20 cm
DAG-206.316	Ø 3 mm X 15 cm	DAG-206.325	Ø 4 mm X 35 cm	DAG-206.333	Ø 6 mm X 25 cm
DAG-206.317	Ø 3 mm X 20 cm	DAG-206.326	Ø 5 mm X 15 cm	DAG-206.334	Ø 6 mm X 30 cm
DAG-206.318	Ø 3 mm X 25 cm	DAG-206.327	Ø 5 mm X 20 cm	DAG-206.335	Ø 6 mm X 35 cm
DAG-206.319	Ø 3 mm X 30 cm				



Bi-Port Cannula


Luer lock (Stainless Steel)

	10cm	15cm	20cm	25cm
2.0 Ø	DAG-205.91	DAG-205.101	DAG-205.111	DAG-205.121
2.5 Ø	DAG-205.92	DAG-205.102	DAG-205.112	DAG-205.122
3.0 Ø	DAG-205.93	DAG-205.103	DAG-205.113	DAG-205.123
3.5 Ø	DAG-205.94	DAG-205.104	DAG-205.114	DAG-205.124
4.0 Ø	DAG-205.95	DAG-205.105	DAG-205.115	DAG-205.125

Power Cannula

DAG-206.351	Ø 2 mm X 15 cm	DAG-206.360	Ø 3 mm X 35 cm	DAG-206.368	Ø 5 mm X 25 cm
DAG-206.352	Ø 2 mm X 20 cm	DAG-206.361	Ø 4 mm X 15 cm	DAG-206.369	Ø 5 mm X 30 cm
DAG-206.353	Ø 2 mm X 25 cm	DAG-206.362	Ø 4 mm X 20 cm	DAG-206.370	Ø 5 mm X 35 cm
DAG-206.354	Ø 2 mm X 30 cm	DAG-206.363	Ø 4 mm X 25 cm	DAG-206.371	Ø 6 mm X 15 cm
DAG-206.355	Ø 2 mm X 35 cm	DAG-206.364	Ø 4 mm X 30 cm	DAG-206.372	Ø 6 mm X 20 cm
DAG-206.356	Ø 3 mm X 15 cm	DAG-206.365	Ø 4 mm X 35 cm	DAG-206.373	Ø 6 mm X 25 cm
DAG-206.357	Ø 3 mm X 20 cm	DAG-206.366	Ø 5 mm X 15 cm	DAG-206.374	Ø 6 mm X 30 cm
DAG-206.358	Ø 3 mm X 25 cm	DAG-206.367	Ø 5 mm X 20 cm	DAG-206.375	Ø 6 mm X 35 cm
DAG-206.359	Ø 3 mm X 30 cm				

General Suction Cannulas



Tappered Tip Single Central Hole

Power Cannula

DAG-206.391	Ø 2 mm X 15 cm	DAG-206.400	Ø 3 mm X 35 cm	DAG-206.408	Ø 5 mm X 25 cm
DAG-206.392	Ø 2 mm X 20 cm	DAG-206.401	Ø 4 mm X 15 cm	DAG-206.409	Ø 5 mm X 30 cm
DAG-206.393	Ø 2 mm X 25 cm	DAG-206.402	Ø 4 mm X 20 cm	DAG-206.410	Ø 5 mm X 35 cm
DAG-206.394	Ø 2 mm X 30 cm	DAG-206.403	Ø 4 mm X 25 cm	DAG-206.411	Ø 6 mm X 15 cm
DAG-206.395	Ø 2 mm X 35 cm	DAG-206.404	Ø 4 mm X 30 cm	DAG-206.412	Ø 6 mm X 20 cm
DAG-206.396	Ø 3 mm X 15 cm	DAG-206.405	Ø 4 mm X 35 cm	DAG-206.413	Ø 6 mm X 25 cm
DAG-206.397	Ø 3 mm X 20 cm	DAG-206.406	Ø 5 mm X 15 cm	DAG-206.414	Ø 6 mm X 30 cm
DAG-206.398	Ø 3 mm X 25 cm	DAG-206.407	Ø 5 mm X 20 cm	DAG-206.415	Ø 6 mm X 35 cm
DAG-206.399	Ø 3 mm X 30 cm				



Tapperd Tip Two Spiral Holes

Power Cannula

DAG-206.431	Ø 2 mm X 15 cm	DAG-206.440	Ø 3 mm X 35 cm	DAG-206.448	Ø 5 mm X 25 cm
DAG-206.432	Ø 2 mm X 20 cm	DAG-206.441	Ø 4 mm X 15 cm	DAG-206.449	Ø 5 mm X 30 cm
DAG-206.433	Ø 2 mm X 25 cm	DAG-206.442	Ø 4 mm X 20 cm	DAG-206.450	Ø 5 mm X 35 cm
DAG-206.434	Ø 2 mm X 30 cm	DAG-206.443	Ø 4 mm X 25 cm	DAG-206.451	Ø 6 mm X 15 cm
DAG-206.435	Ø 2 mm X 35 cm	DAG-206.444	Ø 4 mm X 30 cm	DAG-206.452	Ø 6 mm X 20 cm
DAG-206.436	Ø 3 mm X 15 cm	DAG-206.445	Ø 4 mm X 35 cm	DAG-206.453	Ø 6 mm X 25 cm
DAG-206.437	Ø 3 mm X 20 cm	DAG-206.446	Ø 5 mm X 15 cm	DAG-206.454	Ø 6 mm X 30 cm
DAG-206.438	Ø 3 mm X 25 cm	DAG-206.447	Ø 5 mm X 20 cm	DAG-206.455	Ø 6 mm X 35 cm
DAG-206.439	Ø 3 mm X 30 cm				



Fig. 3F

Power Cannula

DAG-206.471	Ø 2 mm X 15 cm	DAG-206.480	Ø 3 mm X 35 cm	DAG-206.488	Ø 5 mm X 25 cm
DAG-206.472	Ø 2 mm X 20 cm	DAG-206.481	Ø 4 mm X 15 cm	DAG-206.489	Ø 5 mm X 30 cm
DAG-206.473	Ø 2 mm X 25 cm	DAG-206.482	Ø 4 mm X 20 cm	DAG-206.490	Ø 5 mm X 35 cm
DAG-206.474	Ø 2 mm X 30 cm	DAG-206.483	Ø 4 mm X 25 cm	DAG-206.491	Ø 6 mm X 15 cm
DAG-206.475	Ø 2 mm X 35 cm	DAG-206.484	Ø 4 mm X 30 cm	DAG-206.492	Ø 6 mm X 20 cm
DAG-206.476	Ø 3 mm X 15 cm	DAG-206.485	Ø 4 mm X 35 cm	DAG-206.493	Ø 6 mm X 25 cm
DAG-206.477	Ø 3 mm X 20 cm	DAG-206.486	Ø 5 mm X 15 cm	DAG-206.494	Ø 6 mm X 30 cm
DAG-206.478	Ø 3 mm X 25 cm	DAG-206.487	Ø 5 mm X 20 cm	DAG-206.495	Ø 6 mm X 35 cm
DAG-206.479	Ø 3 mm X 30 cm				

General Suction Cannulas




Fig.2Ech

Power Cannula

DAG-206.501	Ø 2 mm X 15 cm	DAG-206.510	Ø 3 mm X 35 cm	DAG-206.518	Ø 5 mm X 25 cm
DAG-206.502	Ø 2 mm X 20 cm	DAG-206.511	Ø 4 mm X 15 cm	DAG-206.519	Ø 5 mm X 30 cm
DAG-206.503	Ø 2 mm X 25 cm	DAG-206.512	Ø 4 mm X 20 cm	DAG-206.520	Ø 5 mm X 35 cm
DAG-206.504	Ø 2 mm X 30 cm	DAG-206.513	Ø 4 mm X 25 cm	DAG-206.521	Ø 6 mm X 15 cm
DAG-206.505	Ø 2 mm X 35 cm	DAG-206.514	Ø 4 mm X 30 cm	DAG-206.522	Ø 6 mm X 20 cm
DAG-206.506	Ø 3 mm X 15 cm	DAG-206.515	Ø 4 mm X 35 cm	DAG-206.523	Ø 6 mm X 25 cm
DAG-206.507	Ø 3 mm X 20 cm	DAG-206.516	Ø 5 mm X 15 cm	DAG-206.524	Ø 6 mm X 30 cm
DAG-206.508	Ø 3 mm X 25 cm	DAG-206.517	Ø 5 mm X 20 cm	DAG-206.525	Ø 6 mm X 35 cm
DAG-206.509	Ø 3 mm X 30 cm				

Two Square-In Line


Luer lock (Stainless Steel)

	7cm	10cm	15cm	20cm	25cm
2.0	ØDAG-205.131	DAG-205.141	DAG-205.151	DAG-205.161	DAG-205.171
2.5	ØDAG-205.132	DAG-205.142	DAG-205.152	DAG-205.162	DAG-205.172
3.0	ØDAG-205.133	DAG-205.143	DAG-205.153	DAG-205.163	DAG-205.173
3.5	ØDAG-205.134	DAG-205.144	DAG-205.154	DAG-205.164	DAG-205.174
4.0	ØDAG-205.135	DAG-205.145	DAG-205.155	DAG-205.165	DAG-205.175

Power Cannula

DAG-206.551	Ø 2 mm X 15 cm	DAG-206.560	Ø 3 mm X 35 cm	DAG-206.568	Ø 5 mm X 25 cm
DAG-206.552	Ø 2 mm X 20 cm	DAG-206.561	Ø 4 mm X 15 cm	DAG-206.569	Ø 5 mm X 30 cm
DAG-206.553	Ø 2 mm X 25 cm	DAG-206.562	Ø 4 mm X 20 cm	DAG-206.570	Ø 5 mm X 35 cm
DAG-206.554	Ø 2 mm X 30 cm	DAG-206.563	Ø 4 mm X 25 cm	DAG-206.571	Ø 6 mm X 15 cm
DAG-206.555	Ø 2 mm X 35 cm	DAG-206.564	Ø 4 mm X 30 cm	DAG-206.572	Ø 6 mm X 20 cm
DAG-206.556	Ø 3 mm X 15 cm	DAG-206.565	Ø 4 mm X 35 cm	DAG-206.573	Ø 6 mm X 25 cm
DAG-206.557	Ø 3 mm X 20 cm	DAG-206.566	Ø 5 mm X 15 cm	DAG-206.574	Ø 6 mm X 30 cm
DAG-206.558	Ø 3 mm X 25 cm	DAG-206.567	Ø 5 mm X 20 cm	DAG-206.575	Ø 6 mm X 35 cm
DAG-206.559	Ø 3 mm X 30 cm				

General Suction Cannulas



One Central Port-Two Parallel Ports

Power Cannula

DAG-206.591	Ø 2 mm X 15 cm	DAG-206.600	Ø 3 mm X 35 cm	DAG-206.608	Ø 5 mm X 25 cm
DAG-206.592	Ø 2 mm X 20 cm	DAG-206.601	Ø 4 mm X 15 cm	DAG-206.609	Ø 5 mm X 30 cm
DAG-206.593	Ø 2 mm X 25 cm	DAG-206.602	Ø 4 mm X 20 cm	DAG-206.610	Ø 5 mm X 35 cm
DAG-206.594	Ø 2 mm X 30 cm	DAG-206.603	Ø 4 mm X 25 cm	DAG-206.611	Ø 6 mm X 15 cm
DAG-206.595	Ø 2 mm X 35 cm	DAG-206.604	Ø 4 mm X 30 cm	DAG-206.612	Ø 6 mm X 20 cm
DAG-206.596	Ø 3 mm X 15 cm	DAG-206.605	Ø 4 mm X 35 cm	DAG-206.613	Ø 6 mm X 25 cm
DAG-206.597	Ø 3 mm X 20 cm	DAG-206.606	Ø 5 mm X 15 cm	DAG-206.614	Ø 6 mm X 30 cm
DAG-206.598	Ø 3 mm X 25 cm	DAG-206.607	Ø 5 mm X 20 cm	DAG-206.615	Ø 6 mm X 35 cm
DAG-206.599	Ø 3 mm X 30 cm				

Two Parallel Port- One Central Port

Power Cannula

DAG-206.651	Ø 2 mm X 15 cm	DAG-206.660	Ø 3 mm X 35 cm	DAG-206.668	Ø 5 mm X 25 cm
DAG-206.652	Ø 2 mm X 20 cm	DAG-206.661	Ø 4 mm X 15 cm	DAG-206.669	Ø 5 mm X 30 cm
DAG-206.653	Ø 2 mm X 25 cm	DAG-206.662	Ø 4 mm X 20 cm	DAG-206.670	Ø 5 mm X 35 cm
DAG-206.654	Ø 2 mm X 30 cm	DAG-206.663	Ø 4 mm X 25 cm	DAG-206.671	Ø 6 mm X 15 cm
DAG-206.655	Ø 2 mm X 35 cm	DAG-206.664	Ø 4 mm X 30 cm	DAG-206.672	Ø 6 mm X 20 cm
DAG-206.656	Ø 3 mm X 15 cm	DAG-206.665	Ø 4 mm X 35 cm	DAG-206.673	Ø 6 mm X 25 cm
DAG-206.657	Ø 3 mm X 20 cm	DAG-206.666	Ø 5 mm X 15 cm	DAG-206.674	Ø 6 mm X 30 cm
DAG-206.658	Ø 3 mm X 25 cm	DAG-206.667	Ø 5 mm X 20 cm	DAG-206.675	Ø 6 mm X 35 cm
DAG-206.659	Ø 3 mm X 30 cm				

General Suction Cannulas




Fig. 3Me

Luer lock (Stainless Steel)

	7cm	10cm	15cm	20cm	25cm
2.0	Ø DAG-205.181	DAG-205.191	DAG-205.201	DAG-205.211	DAG-205.221
2.5	Ø DAG-205.182	DAG-205.192	DAG-205.202	DAG-205.212	DAG-205.222
3.0	Ø DAG-205.183	DAG-205.193	DAG-205.203	DAG-205.213	DAG-205.223
3.5	Ø DAG-205.184	DAG-205.194	DAG-205.204	DAG-205.214	DAG-205.224
4.0	Ø DAG-205.185	DAG-205.195	DAG-205.205	DAG-205.215	DAG-205.225

Power Cannula

DAG-206.681	Ø 2 mm X 15 cm	DAG-206.690	Ø 3 mm X 35 cm	DAG-206.698	Ø 5 mm X 25 cm
DAG-206.682	Ø 2 mm X 20 cm	DAG-206.691	Ø 4 mm X 15 cm	DAG-206.699	Ø 5 mm X 30 cm
DAG-206.683	Ø 2 mm X 25 cm	DAG-206.692	Ø 4 mm X 20 cm	DAG-206.700	Ø 5 mm X 35 cm
DAG-206.684	Ø 2 mm X 30 cm	DAG-206.693	Ø 4 mm X 25 cm	DAG-206.701	Ø 6 mm X 15 cm
DAG-206.685	Ø 2 mm X 35 cm	DAG-206.694	Ø 4 mm X 30 cm	DAG-206.702	Ø 6 mm X 20 cm
DAG-206.686	Ø 3 mm X 15 cm	DAG-206.695	Ø 4 mm X 35 cm	DAG-206.703	Ø 6 mm X 25 cm
DAG-206.687	Ø 3 mm X 20 cm	DAG-206.696	Ø 5 mm X 15 cm	DAG-206.704	Ø 6 mm X 30 cm
DAG-206.688	Ø 3 mm X 25 cm	DAG-206.697	Ø 5 mm X 20 cm	DAG-206.705	Ø 6 mm X 35 cm
DAG-206.689	Ø 3 mm X 30 cm				

Fig. 3Hr

Power Cannula

DAG-206.721	Ø 2 mm X 15 cm	DAG-206.730	Ø 3 mm X 35 cm	DAG-206.738	Ø 5 mm X 25 cm
DAG-206.722	Ø 2 mm X 20 cm	DAG-206.731	Ø 4 mm X 15 cm	DAG-206.739	Ø 5 mm X 30 cm
DAG-206.723	Ø 2 mm X 25 cm	DAG-206.732	Ø 4 mm X 20 cm	DAG-206.740	Ø 5 mm X 35 cm
DAG-206.724	Ø 2 mm X 30 cm	DAG-206.733	Ø 4 mm X 25 cm	DAG-206.741	Ø 6 mm X 15 cm
DAG-206.725	Ø 2 mm X 35 cm	DAG-206.734	Ø 4 mm X 30 cm	DAG-206.742	Ø 6 mm X 20 cm
DAG-206.726	Ø 3 mm X 15 cm	DAG-206.735	Ø 4 mm X 35 cm	DAG-206.743	Ø 6 mm X 25 cm
DAG-206.727	Ø 3 mm X 20 cm	DAG-206.736	Ø 5 mm X 15 cm	DAG-206.744	Ø 6 mm X 30 cm
DAG-206.728	Ø 3 mm X 25 cm	DAG-206.737	Ø 5 mm X 20 cm	DAG-206.745	Ø 6 mm X 35 cm
DAG-206.729	Ø 3 mm X 30 cm				

Fig. Riz

Power Cannula

DAG-206.751	Ø 2 mm X 15 cm	DAG-206.760	Ø 3 mm X 35 cm	DAG-206.768	Ø 5 mm X 25 cm
DAG-206.752	Ø 2 mm X 20 cm	DAG-206.761	Ø 4 mm X 15 cm	DAG-206.769	Ø 5 mm X 30 cm
DAG-206.753	Ø 2 mm X 25 cm	DAG-206.762	Ø 4 mm X 20 cm	DAG-206.770	Ø 5 mm X 35 cm
DAG-206.754	Ø 2 mm X 30 cm	DAG-206.763	Ø 4 mm X 25 cm	DAG-206.771	Ø 6 mm X 15 cm
DAG-206.755	Ø 2 mm X 35 cm	DAG-206.764	Ø 4 mm X 30 cm	DAG-206.772	Ø 6 mm X 20 cm
DAG-206.756	Ø 3 mm X 15 cm	DAG-206.765	Ø 4 mm X 35 cm	DAG-206.773	Ø 6 mm X 25 cm
DAG-206.757	Ø 3 mm X 20 cm	DAG-206.766	Ø 5 mm X 15 cm	DAG-206.774	Ø 6 mm X 30 cm
DAG-206.758	Ø 3 mm X 25 cm	DAG-206.767	Ø 5 mm X 20 cm	DAG-206.775	Ø 6 mm X 35 cm
DAG-206.759	Ø 3 mm X 30 cm				

General Suction Cannulas




Fig. 5Ai

Power Cannula

DAG-206.721	Ø 2 mm X 15 cm	DAG-206.730	Ø 3 mm X 35 cm	DAG-206.738	Ø 5 mm X 25 cm
DAG-206.722	Ø 2 mm X 20 cm	DAG-206.731	Ø 4 mm X 15 cm	DAG-206.739	Ø 5 mm X 30 cm
DAG-206.723	Ø 2 mm X 25 cm	DAG-206.732	Ø 4 mm X 20 cm	DAG-206.740	Ø 5 mm X 35 cm
DAG-206.724	Ø 2 mm X 30 cm	DAG-206.733	Ø 4 mm X 25 cm	DAG-206.741	Ø 6 mm X 15 cm
DAG-206.725	Ø 2 mm X 35 cm	DAG-206.734	Ø 4 mm X 30 cm	DAG-206.742	Ø 6 mm X 20 cm
DAG-206.726	Ø 3 mm X 15 cm	DAG-206.735	Ø 4 mm X 35 cm	DAG-206.743	Ø 6 mm X 25 cm
DAG-206.727	Ø 3 mm X 20 cm	DAG-206.736	Ø 5 mm X 15 cm	DAG-206.744	Ø 6 mm X 30 cm
DAG-206.728	Ø 3 mm X 25 cm	DAG-206.737	Ø 5 mm X 20 cm	DAG-206.745	Ø 6 mm X 35 cm
DAG-206.729	Ø 3 mm X 30 cm				

Fig. Mr

Power Cannula

DAG-206.751	Ø 2 mm X 15 cm	DAG-206.760	Ø 3 mm X 35 cm	DAG-206.768	Ø 5 mm X 25 cm
DAG-206.752	Ø 2 mm X 20 cm	DAG-206.761	Ø 4 mm X 15 cm	DAG-206.769	Ø 5 mm X 30 cm
DAG-206.753	Ø 2 mm X 25 cm	DAG-206.762	Ø 4 mm X 20 cm	DAG-206.770	Ø 5 mm X 35 cm
DAG-206.754	Ø 2 mm X 30 cm	DAG-206.763	Ø 4 mm X 25 cm	DAG-206.771	Ø 6 mm X 15 cm
DAG-206.755	Ø 2 mm X 35 cm	DAG-206.764	Ø 4 mm X 30 cm	DAG-206.772	Ø 6 mm X 20 cm
DAG-206.756	Ø 3 mm X 15 cm	DAG-206.765	Ø 4 mm X 35 cm	DAG-206.773	Ø 6 mm X 25 cm
DAG-206.757	Ø 3 mm X 20 cm	DAG-206.766	Ø 5 mm X 15 cm	DAG-206.774	Ø 6 mm X 30 cm
DAG-206.758	Ø 3 mm X 25 cm	DAG-206.767	Ø 5 mm X 20 cm	DAG-206.775	Ø 6 mm X 35 cm
DAG-206.759	Ø 3 mm X 30 cm				

Fig. P3G

Power Cannula

DAG-206.791	Ø 2 mm X 15 cm	DAG-206.800	Ø 3 mm X 35 cm	DAG-206.808	Ø 5 mm X 25 cm
DAG-206.792	Ø 2 mm X 20 cm	DAG-206.801	Ø 4 mm X 15 cm	DAG-206.809	Ø 5 mm X 30 cm
DAG-206.793	Ø 2 mm X 25 cm	DAG-206.802	Ø 4 mm X 20 cm	DAG-206.810	Ø 5 mm X 35 cm
DAG-206.794	Ø 2 mm X 30 cm	DAG-206.803	Ø 4 mm X 25 cm	DAG-206.811	Ø 6 mm X 15 cm
DAG-206.795	Ø 2 mm X 35 cm	DAG-206.804	Ø 4 mm X 30 cm	DAG-206.812	Ø 6 mm X 20 cm
DAG-206.796	Ø 3 mm X 15 cm	DAG-206.805	Ø 4 mm X 35 cm	DAG-206.813	Ø 6 mm X 25 cm
DAG-206.797	Ø 3 mm X 20 cm	DAG-206.806	Ø 5 mm X 15 cm	DAG-206.814	Ø 6 mm X 30 cm
DAG-206.798	Ø 3 mm X 25 cm	DAG-206.807	Ø 5 mm X 20 cm	DAG-206.815	Ø 6 mm X 35 cm
DAG-206.799	Ø 3 mm X 30 cm				

General Suction Cannulas




Fig.Py

Power Cannula

DAG-206.831	Ø 2 mm X 15 cm	DAG-206.840	Ø 3 mm X 35 cm	DAG-206.848	Ø 5 mm X 25 cm
DAG-206.832	Ø 2 mm X 20 cm	DAG-206.841	Ø 4 mm X 15 cm	DAG-206.849	Ø 5 mm X 30 cm
DAG-206.833	Ø 2 mm X 25 cm	DAG-206.842	Ø 4 mm X 20 cm	DAG-206.850	Ø 5 mm X 35 cm
DAG-206.834	Ø 2 mm X 30 cm	DAG-206.843	Ø 4 mm X 25 cm	DAG-206.851	Ø 6 mm X 15 cm
DAG-206.835	Ø 2 mm X 35 cm	DAG-206.844	Ø 4 mm X 30 cm	DAG-206.852	Ø 6 mm X 20 cm
DAG-206.836	Ø 3 mm X 15 cm	DAG-206.845	Ø 4 mm X 35 cm	DAG-206.853	Ø 6 mm X 25 cm
DAG-206.837	Ø 3 mm X 20 cm	DAG-206.846	Ø 5 mm X 15 cm	DAG-206.854	Ø 6 mm X 30 cm
DAG-206.838	Ø 3 mm X 25 cm	DAG-206.847	Ø 5 mm X 20 cm	DAG-206.855	Ø 6 mm X 35 cm
DAG-206.839	Ø 3 mm X 30 cm				

Fig.Le


Power Cannula

DAG-206.861	Ø 2 mm X 15 cm	DAG-206.870	Ø 3 mm X 35 cm	DAG-206.878	Ø 5 mm X 25 cm
DAG-206.862	Ø 2 mm X 20 cm	DAG-206.871	Ø 4 mm X 15 cm	DAG-206.879	Ø 5 mm X 30 cm
DAG-206.863	Ø 2 mm X 25 cm	DAG-206.872	Ø 4 mm X 20 cm	DAG-206.880	Ø 5 mm X 35 cm
DAG-206.864	Ø 2 mm X 30 cm	DAG-206.873	Ø 4 mm X 25 cm	DAG-206.881	Ø 6 mm X 15 cm
DAG-206.865	Ø 2 mm X 35 cm	DAG-206.874	Ø 4 mm X 30 cm	DAG-206.882	Ø 6 mm X 20 cm
DAG-206.866	Ø 3 mm X 15 cm	DAG-206.875	Ø 4 mm X 35 cm	DAG-206.883	Ø 6 mm X 25 cm
DAG-206.867	Ø 3 mm X 20 cm	DAG-206.876	Ø 5 mm X 15 cm	DAG-206.884	Ø 6 mm X 30 cm
DAG-206.868	Ø 3 mm X 25 cm	DAG-206.877	Ø 5 mm X 20 cm	DAG-206.885	Ø 6 mm X 35 cm
DAG-206.869	Ø 3 mm X 30 cm				

HARVESTER Cannulas

Fat Harvesting for Re-grafting is a delicate difficult and time consuming process. But with Dag Surgical corp's **Designedly** harvesting Cannulas will help you to achieve this goal in super efficient way.

Micro fat Harvesting - No problem with Dag Surgical corporation Cannulas.



Harvester Cannulas



Multiport Speed Harvester- Both End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

60ccTommeyhub

DAG-301.01	Ø 2.5 mm X 15 cm
DAG-301.02	Ø 2.5 mm X 20 cm
DAG-301.03	Ø 3.0 mm X 15 cm
DAG-301.04	Ø 3.0 mm X 20 cm
DAG-301.05	Ø 3.0 mm X 25 cm
DAG-301.06	Ø 4.0 mm X 25 cm
DAG-301.07	Ø 4.0 mm X 30 cm
DAG-301.08	Ø 4.0 mm X 35 cm
DAG-301.09	Ø 4.5 mm X 30 cm
DAG-301.10	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-302.01	Ø 3.0 mm X 15 cm
DAG-302.02	Ø 3.0 mm X 25 cm
DAG-302.03	Ø 4.0 mm X 15 cm
DAG-302.04	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-303.01	Ø 2.0 mm X 15 cm
DAG-303.02	Ø 2.5 mm X 15 cm
DAG-303.03	Ø 2.5 mm X 20 cm
DAG-303.04	Ø 3.0 mm X 15 cm
DAG-303.05	Ø 3.0 mm X 20 cm
DAG-303.06	Ø 3.0 mm X 25 cm

Luer lock

DAG-304.01	Ø 2.0 mm X 15 cm
DAG-304.02	Ø 2.5 mm X 15 cm
DAG-304.03	Ø 2.5 mm X 20 cm
DAG-304.04	Ø 3.0 mm X 15 cm
DAG-304.05	Ø 3.0 mm X 20 cm
DAG-304.06	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-305.01	Ø 2.0 mm X 15 cm
DAG-305.02	Ø 2.5 mm X 15 cm
DAG-305.03	Ø 2.5 mm X 20 cm
DAG-305.04	Ø 3.0 mm X 15 cm
DAG-305.05	Ø 3.0 mm X 20 cm
DAG-305.06	Ø 3.0 mm X 25 cm

Power Cannula

DAG-306.01	Ø 2 mm X 15 cm
DAG-306.02	Ø 2 mm X 20 cm
DAG-306.03	Ø 2 mm X 25 cm
DAG-306.04	Ø 2 mm X 30 cm
DAG-306.05	Ø 2 mm X 35 cm
DAG-306.06	Ø 3 mm X 15 cm
DAG-306.07	Ø 3 mm X 20 cm
DAG-306.08	Ø 3 mm X 25 cm
DAG-306.09	Ø 3 mm X 30 cm
DAG-306.10	Ø 3 mm X 35 cm
DAG-306.11	Ø 4 mm X 15 cm
DAG-306.12	Ø 4 mm X 20 cm
DAG-306.13	Ø 4 mm X 25 cm
DAG-306.14	Ø 4 mm X 30 cm
DAG-306.15	Ø 4 mm X 35 cm
DAG-306.16	Ø 5 mm X 15 cm
DAG-306.17	Ø 5 mm X 20 cm
DAG-306.18	Ø 5 mm X 25 cm
DAG-306.19	Ø 5 mm X 30 cm
DAG-306.20	Ø 5 mm X 35 cm
DAG-306.21	Ø 6 mm X 15 cm
DAG-306.22	Ø 6 mm X 20 cm
DAG-306.23	Ø 6 mm X 25 cm
DAG-306.24	Ø 6 mm X 30 cm
DAG-306.25	Ø 6 mm X 35 cm

Harvester Cannulas



Multiport Speed Harvester- One End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

60ccTommeyhub

DAG-301.11	Ø 2.5 mm X 15 cm
DAG-301.12	Ø 2.5 mm X 20 cm
DAG-301.13	Ø 3.0 mm X 15 cm
DAG-301.14	Ø 3.0 mm X 20 cm
DAG-301.15	Ø 3.0 mm X 25 cm
DAG-301.16	Ø 4.0 mm X 25 cm
DAG-301.17	Ø 4.0 mm X 30 cm
DAG-301.18	Ø 4.0 mm X 35 cm
DAG-301.19	Ø 4.5 mm X 30 cm
DAG-301.20	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-302.11	Ø 3.0 mm X 15 cm
DAG-302.12	Ø 3.0 mm X 25 cm
DAG-302.13	Ø 4.0 mm X 15 cm
DAG-302.14	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-303.11	Ø 2.0 mm X 15 cm
DAG-303.12	Ø 2.5 mm X 15 cm
DAG-303.13	Ø 2.5 mm X 20 cm
DAG-303.14	Ø 3.0 mm X 15 cm
DAG-303.15	Ø 3.0 mm X 20 cm
DAG-303.16	Ø 3.0 mm X 25 cm

Luer lock

DAG-304.11	Ø 2.0 mm X 15 cm
DAG-304.12	Ø 2.5 mm X 15 cm
DAG-304.13	Ø 2.5 mm X 20 cm
DAG-304.14	Ø 3.0 mm X 15 cm
DAG-304.15	Ø 3.0 mm X 20 cm
DAG-304.16	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-305.11	Ø 2.0 mm X 15 cm
DAG-305.12	Ø 2.5 mm X 15 cm
DAG-305.13	Ø 2.5 mm X 20 cm
DAG-305.14	Ø 3.0 mm X 15 cm
DAG-305.15	Ø 3.0 mm X 20 cm
DAG-305.16	Ø 3.0 mm X 25 cm

Power Cannula

DAG-306.11	Ø 2 mm X 15 cm
DAG-306.12	Ø 2 mm X 20 cm
DAG-306.13	Ø 2 mm X 25 cm
DAG-306.14	Ø 2 mm X 30 cm
DAG-306.15	Ø 2 mm X 35 cm
DAG-306.16	Ø 3 mm X 15 cm
DAG-306.17	Ø 3 mm X 20 cm
DAG-306.18	Ø 3 mm X 25 cm
DAG-306.19	Ø 3 mm X 30 cm
DAG-306.20	Ø 3 mm X 35 cm
DAG-306.21	Ø 4 mm X 15 cm
DAG-306.22	Ø 4 mm X 20 cm
DAG-306.23	Ø 4 mm X 25 cm
DAG-306.24	Ø 4 mm X 30 cm
DAG-306.25	Ø 4 mm X 35 cm
DAG-306.26	Ø 5 mm X 15 cm
DAG-306.27	Ø 5 mm X 20 cm
DAG-306.28	Ø 5 mm X 25 cm
DAG-306.29	Ø 5 mm X 30 cm
DAG-306.30	Ø 5 mm X 35 cm
DAG-306.31	Ø 6 mm X 15 cm
DAG-306.32	Ø 6 mm X 20 cm
DAG-306.33	Ø 6 mm X 25 cm
DAG-306.34	Ø 6 mm X 30 cm
DAG-306.35	Ø 6 mm X 35 cm

Harvester Cannulas



Multi Port Harvester - Both End Bevel

A sophisticated, precision-designed harvester that is truly Patient Friendly. The ports in the Sorensen Harvester are machined with alternating 60 degree cutting edges. This unique design is less traumatic than other multi-port harvesters and delivers a fine-grained, high fat content/purity yield that can easily pass through 1.2mm, 0.9mm and 0.7mm blunt injectors.

60ccTommeYhub

DAG-301.21	Ø 2.5 mm X 15 cm
DAG-301.22	Ø 2.5 mm X 20 cm
DAG-301.23	Ø 3.0 mm X 15 cm
DAG-301.24	Ø 3.0 mm X 20 cm
DAG-301.25	Ø 3.0 mm X 25 cm
DAG-301.26	Ø 4.0 mm X 25 cm
DAG-301.27	Ø 4.0 mm X 30 cm
DAG-301.28	Ø 4.0 mm X 35 cm
DAG-301.29	Ø 4.5 mm X 30 cm
DAG-301.30	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-302.21	Ø 3.0 mm X 15 cm
DAG-302.22	Ø 3.0 mm X 25 cm
DAG-302.23	Ø 4.0 mm X 15 cm
DAG-302.24	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-303.21	Ø 2.0 mm X 15 cm
DAG-303.22	Ø 2.5 mm X 15 cm
DAG-303.23	Ø 2.5 mm X 20 cm
DAG-303.24	Ø 3.0 mm X 15 cm
DAG-303.25	Ø 3.0 mm X 20 cm
DAG-303.26	Ø 3.0 mm X 25 cm

Luer lock

DAG-304.21	Ø 2.0 mm X 15 cm
DAG-304.22	Ø 2.5 mm X 15 cm
DAG-304.23	Ø 2.5 mm X 20 cm
DAG-304.24	Ø 3.0 mm X 15 cm
DAG-304.25	Ø 3.0 mm X 20 cm
DAG-304.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)


DAG-305.21	Ø 2.0 mm X 15 cm
DAG-305.22	Ø 2.5 mm X 15 cm
DAG-305.23	Ø 2.5 mm X 20 cm
DAG-305.24	Ø 3.0 mm X 15 cm
DAG-305.25	Ø 3.0 mm X 20 cm
DAG-305.26	Ø 3.0 mm X 25 cm

Power Cannula

DAG-306.21	Ø 2 mm X 15 cm
DAG-306.22	Ø 2 mm X 20 cm
DAG-306.23	Ø 2 mm X 25 cm
DAG-306.24	Ø 2 mm X 30 cm
DAG-306.25	Ø 2 mm X 35 cm
DAG-306.26	Ø 3 mm X 15 cm
DAG-306.27	Ø 3 mm X 20 cm
DAG-306.28	Ø 3 mm X 25 cm
DAG-306.29	Ø 3 mm X 30 cm
DAG-306.30	Ø 3 mm X 35 cm
DAG-306.31	Ø 4 mm X 15 cm
DAG-306.32	Ø 4 mm X 20 cm
DAG-306.33	Ø 4 mm X 25 cm
DAG-306.34	Ø 4 mm X 30 cm
DAG-306.35	Ø 4 mm X 35 cm
DAG-306.36	Ø 5 mm X 15 cm
DAG-306.37	Ø 5 mm X 20 cm
DAG-306.38	Ø 5 mm X 25 cm
DAG-306.39	Ø 5 mm X 30 cm
DAG-306.40	Ø 5 mm X 35 cm
DAG-306.41	Ø 6 mm X 15 cm
DAG-306.42	Ø 6 mm X 20 cm
DAG-306.43	Ø 6 mm X 25 cm
DAG-306.44	Ø 6 mm X 30 cm
DAG-306.45	Ø 6 mm X 35 cm

INJECTOR Cannulas

(Micro Injectors)



Injector Cannulas



Sharp Extractor/Injector

This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

Ultra Luer lock		Luer lock	
DAG-403.01	Ø 1.0 mm X 5 cm	DAG-404.01	Ø 1.0 mm X 5 cm
DAG-403.02	Ø 1.0 mm X 7 cm	DAG-404.02	Ø 1.0 mm X 7 cm
DAG-403.03	Ø 2.0 mm X 15 cm	DAG-404.03	Ø 2.0 mm X 15 cm
DAG-403.04	Ø 2.5 mm X 15 cm	DAG-404.04	Ø 2.5 mm X 15 cm
DAG-403.05	Ø 3.0 mm X 15 cm	DAG-404.05	Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.01	DAG-405.11	DAG-405.21	DAG-405.31	DAG-405.41
1.4Ø	DAG-405.02	DAG-405.12	DAG-405.22	DAG-405.32	DAG-405.42
1.8Ø	DAG-405.03	DAG-405.13	DAG-405.23	DAG-405.33	DAG-405.43
2.0Ø	DAG-405.04	DAG-405.14	DAG-405.24	DAG-405.34	DAG-405.44

Blunt Extractor/Injector

This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

Ultra Luer lock		Luer lock	
DAG-403.11	Ø 2.0 mm X 15 cm	DAG-404.11	Ø 2.0 mm X 15 cm
DAG-403.12	Ø 2.5 mm X 15 cm	DAG-404.12	Ø 2.5 mm X 15 cm
DAG-403.13	Ø 3.0 mm X 15 cm	DAG-404.13	Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.51	DAG-405.61	DAG-405.71	DAG-405.81	DAG-405.91
1.4Ø	DAG-405.52	DAG-405.62	DAG-405.72	DAG-405.82	DAG-405.92
1.8Ø	DAG-405.53	DAG-405.63	DAG-405.73	DAG-405.83	DAG-405.93
2.0Ø	DAG-405.54	DAG-405.64	DAG-405.74	DAG-405.84	DAG-405.94

Blunt Tip Injector

Rounded blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp instruments.

Ultra Luer lock		Luer lock	
DAG-403.21	Ø 1.0 mm X 4 cm	DAG-404.21	Ø 1.0 mm X 4 cm
DAG-403.22	Ø 1.0 mm X 5 cm	DAG-404.22	Ø 1.0 mm X 5 cm
DAG-403.23	Ø 1.0 mm X 7 cm	DAG-404.23	Ø 1.0 mm X 7 cm
DAG-403.24	Curved 12 cm	DAG-404.24	Curved 12 cm
DAG-403.25	Curved 15 cm	DAG-404.25	Curved 15 cm
DAG-403.26	Curved 20 cm	DAG-404.26	Curved 20 cm

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.101	DAG-405.111	DAG-405.121	DAG-405.131	DAG-405.141
1.4Ø	DAG-405.102	DAG-405.112	DAG-405.122	DAG-405.132	DAG-405.142
1.8Ø	DAG-405.103	DAG-405.113	DAG-405.123	DAG-405.133	DAG-405.143
2.0Ø	DAG-405.104	DAG-405.114	DAG-405.124	DAG-405.134	DAG-405.144

Injector Cannulas




Fig. 2B

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.201	DAG-405.211	DAG-405.221	DAG-405.231	DAG-405.241
1.4Ø	DAG-405.202	DAG-405.212	DAG-405.222	DAG-405.232	DAG-405.242
1.8Ø	DAG-405.203	DAG-405.213	DAG-405.223	DAG-405.233	DAG-405.243
2.0Ø	DAG-405.204	DAG-405.214	DAG-405.224	DAG-405.234	DAG-405.244

Fig. 1C

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.251	DAG-405.261	DAG-405.271	DAG-405.281	DAG-405.291
1.4Ø	DAG-405.252	DAG-405.262	DAG-405.272	DAG-405.282	DAG-405.292
1.8Ø	DAG-405.253	DAG-405.263	DAG-405.273	DAG-405.283	DAG-405.293
2.0Ø	DAG-405.254	DAG-405.264	DAG-405.274	DAG-405.284	DAG-405.294


Fig. 1P

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.301	DAG-405.311	DAG-405.321	DAG-405.331	DAG-405.341
1.4Ø	DAG-405.302	DAG-405.312	DAG-405.322	DAG-405.332	DAG-405.342
1.8Ø	DAG-405.303	DAG-405.313	DAG-405.323	DAG-405.333	DAG-405.343
2.0Ø	DAG-405.304	DAG-405.314	DAG-405.324	DAG-405.334	DAG-405.344

Fig. 2E

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.351	DAG-405.361	DAG-405.371	DAG-405.381	DAG-405.391
1.4Ø	DAG-405.352	DAG-405.362	DAG-405.372	DAG-405.382	DAG-405.392
1.8Ø	DAG-405.353	DAG-405.363	DAG-405.373	DAG-405.383	DAG-405.393
2.0Ø	DAG-405.354	DAG-405.364	DAG-405.374	DAG-405.384	DAG-405.394

Injector Cannulas



 Coleman Liposuction Cannula Fig. MIC I

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm	3.5cm
1.1Ø	DAG-405.401	DAG-405.411	DAG-405.421	DAG-405.431	DAG-405.441	DAG-405.451	DAG-405.461	DAG-405.471
1.4Ø	DAG-405.402	DAG-405.412	DAG-405.422	DAG-405.432	DAG-405.442	DAG-405.452	DAG-405.462	DAG-405.472
1.8Ø	DAG-405.403	DAG-405.413	DAG-405.423	DAG-405.433	DAG-405.443	DAG-405.453	DAG-405.463	DAG-405.473
2.0Ø	DAG-405.404	DAG-405.414	DAG-405.424	DAG-405.434	DAG-405.444	DAG-405.454	DAG-405.464	DAG-405.474
2.5Ø	DAG-405.405	DAG-405.415	DAG-405.425	DAG-405.435	DAG-405.445	DAG-405.455	DAG-405.465	DAG-405.475

 Coleman Liposuction Cannula Fig. MIC II

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	DAG-405.501	DAG-405.511	DAG-405.521	DAG-405.531	DAG-405.541	DAG-405.551	DAG-405.561
1.4Ø	DAG-405.502	DAG-405.512	DAG-405.522	DAG-405.532	DAG-405.542	DAG-405.552	DAG-405.562
1.8Ø	DAG-405.503	DAG-405.513	DAG-405.523	DAG-405.533	DAG-405.543	DAG-405.553	DAG-405.563
2.0Ø	DAG-405.504	DAG-405.514	DAG-405.524	DAG-405.534	DAG-405.544	DAG-405.554	DAG-405.564
2.5Ø	DAG-405.505	DAG-405.515	DAG-405.525	DAG-405.535	DAG-405.545	DAG-405.555	DAG-405.565

 Coleman Liposuction Cannula Spoon Tip Fig. MIC III

This tip works well for needle extraction of soft tissue and as a soft tissue re injection needle.

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	DAG-405.601	DAG-405.611	DAG-405.621	DAG-405.631	DAG-405.641	DAG-405.651	DAG-405.661
1.4Ø	DAG-405.602	DAG-405.612	DAG-405.622	DAG-405.632	DAG-405.642	DAG-405.652	DAG-405.662
1.8Ø	DAG-405.603	DAG-405.613	DAG-405.623	DAG-405.633	DAG-405.643	DAG-405.653	DAG-405.663
2.0Ø	DAG-405.604	DAG-405.614	DAG-405.624	DAG-405.634	DAG-405.644	DAG-405.654	DAG-405.664
2.5Ø	DAG-405.605	DAG-405.615	DAG-405.625	DAG-405.635	DAG-405.645	DAG-405.655	DAG-405.665

Little Injector

 Open and more pointed blunt tip is very effective in dense and fibrous tissue.

Ultra Luer lock

DAG-403.51	Ø 1.0 mm X 4 cm
DAG-403.52	Ø 1.0 mm X 5 cm
DAG-403.53	Ø 1.0 mm X 7 cm

Luer lock

DAG-404.51	Ø 1.0 mm X 4 cm
DAG-404.52	Ø 1.0 mm X 5 cm
DAG-404.53	Ø 1.0 mm X 7 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	DAG-405.701	DAG-405.711	DAG-405.721	DAG-405.731	DAG-405.741
1.4Ø	DAG-405.702	DAG-405.712	DAG-405.722	DAG-405.732	DAG-405.742
1.8Ø	DAG-405.703	DAG-405.713	DAG-405.723	DAG-405.733	DAG-405.743
2.0Ø	DAG-405.704	DAG-405.714	DAG-405.724	DAG-405.734	DAG-405.744

Injector Cannulas




 Fig. 1B

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	DAG-405.801	DAG-405.811	DAG-405.821	DAG-405.831	DAG-405.841	DAG-405.851	DAG-405.861
1.4Ø	DAG-405.802	DAG-405.812	DAG-405.822	DAG-405.832	DAG-405.842	DAG-405.852	DAG-405.862
1.8Ø	DAG-405.803	DAG-405.813	DAG-405.823	DAG-405.833	DAG-405.843	DAG-405.853	DAG-405.863
2.0Ø	DAG-405.804	DAG-405.814	DAG-405.824	DAG-405.834	DAG-405.844	DAG-405.854	DAG-405.864
2.5Ø	DAG-405.805	DAG-405.815	DAG-405.825	DAG-405.835	DAG-405.845	DAG-405.855	DAG-405.865

Curved Injectors



Curved Injectors

New Curved Injectors are designed especially for fat transfer procedures to the body. Their unique shape mimics the anatomical curvature of the female torso and, thus, easily produces the natural-looking results that patients desire. Curved Injection Cannulas ameliorate the fat transfer process by providing a direct and precise route to the injection site and by lessening the time required to correctly sculpt curves.

SPECIALITY Cannulas

DAG SURGICAL CORPORATION Speciality Cannula are useful tools for different to operation are such as **face** or **Neck**.

Body

SHAPE YOUR

with




DAG

SURGICAL CORPORATION

4647



Speciality Cannulas



V-Dissector

This cannula is specifically designed for use in the superficial compartment of the truncal area. The cutting notch is very useful for undermining in facial procedures and in dense, fibrous tissue.

60ccTommeyhub	
DAG-501.01	Ø 2.5 mm X 15 cm
DAG-501.02	Ø 2.5 mm X 20 cm
DAG-501.03	Ø 3.0 mm X 15 cm
DAG-501.04	Ø 3.0 mm X 20 cm
DAG-501.05	Ø 3.0 mm X 25 cm
DAG-501.06	Ø 4.0 mm X 25 cm
DAG-501.07	Ø 4.0 mm X 30 cm
DAG-501.08	Ø 4.0 mm X 35 cm
DAG-501.09	Ø 4.5 mm X 30 cm
DAG-501.10	Ø 4.5 mm X 35 cm

Treaded Hub	
DAG-502.01	Ø 3.0 mm X 15 cm
DAG-502.02	Ø 3.0 mm X 25 cm
DAG-502.03	Ø 4.0 mm X 15 cm
DAG-502.04	Ø 4.0 mm X 25 cm


Ultra Luer lock	
DAG-503.01	Ø 2.0 mm X 15 cm
DAG-503.02	Ø 2.5 mm X 15 cm
DAG-503.03	Ø 2.5 mm X 20 cm
DAG-503.04	Ø 3.0 mm X 15 cm
DAG-503.05	Ø 3.0 mm X 20 cm
DAG-503.06	Ø 3.0 mm X 25 cm

Luer lock	
DAG-504.01	Ø 2.0 mm X 15 cm
DAG-504.02	Ø 2.5 mm X 15 cm
DAG-504.03	Ø 2.5 mm X 20 cm
DAG-504.04	Ø 3.0 mm X 15 cm
DAG-504.05	Ø 3.0 mm X 20 cm
DAG-504.06	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)	
DAG-505.01	Ø 2.0 mm X 15 cm
DAG-505.02	Ø 2.5 mm X 15 cm
DAG-505.03	Ø 2.5 mm X 20 cm
DAG-505.04	Ø 3.0 mm X 15 cm
DAG-505.05	Ø 3.0 mm X 20 cm
DAG-505.06	Ø 3.0 mm X 25 cm

Power Cannula	
DAG-506.01	Ø 2 mm X 15 cm
DAG-506.02	Ø 2 mm X 20 cm
DAG-506.03	Ø 2 mm X 25 cm
DAG-506.04	Ø 2 mm X 30 cm
DAG-506.05	Ø 2 mm X 35 cm
DAG-506.06	Ø 3 mm X 15 cm
DAG-506.07	Ø 3 mm X 20 cm
DAG-506.08	Ø 3 mm X 25 cm
DAG-506.09	Ø 3 mm X 30 cm
DAG-506.10	Ø 3 mm X 35 cm
DAG-506.11	Ø 4 mm X 15 cm
DAG-506.12	Ø 4 mm X 20 cm
DAG-506.13	Ø 4 mm X 25 cm
DAG-506.14	Ø 4 mm X 30 cm
DAG-506.15	Ø 4 mm X 35 cm
DAG-506.16	Ø 5 mm X 15 cm
DAG-506.17	Ø 5 mm X 20 cm
DAG-506.18	Ø 5 mm X 25 cm
DAG-506.19	Ø 5 mm X 30 cm
DAG-506.20	Ø 5 mm X 35 cm
DAG-506.21	Ø 6 mm X 15 cm
DAG-506.22	Ø 6 mm X 20 cm
DAG-506.23	Ø 6 mm X 25 cm
DAG-506.24	Ø 6 mm X 30 cm
DAG-506.25	Ø 6 mm X 35 cm

Speciality Cannulas



Flap

The Flap is primarily used as an undermining dissector, and in smaller diameters, for facial tissue removal. The Flap Cannula in a 2.1mm diameter is very useful for undermining the nasolabial fold prior to autologous tissue injections, as well as for tissue removal in this anatomic area.

60ccTommeyhub	
DAG-501.11	Ø 2.5 mm X 15 cm
DAG-501.12	Ø 2.5 mm X 20 cm
DAG-501.13	Ø 3.0 mm X 15 cm
DAG-501.14	Ø 3.0 mm X 20 cm
DAG-501.15	Ø 3.0 mm X 25 cm
DAG-501.16	Ø 4.0 mm X 25 cm
DAG-501.17	Ø 4.0 mm X 30 cm
DAG-501.18	Ø 4.0 mm X 35 cm
DAG-501.19	Ø 4.5 mm X 30 cm
DAG-501.20	Ø 4.5 mm X 35 cm

Treaded Hub	
DAG-502.11	Ø 3.0 mm X 15 cm
DAG-502.12	Ø 3.0 mm X 25 cm
DAG-502.13	Ø 4.0 mm X 15 cm
DAG-502.14	Ø 4.0 mm X 25 cm


Ultra Luer lock	
DAG-503.11	Ø 2.0 mm X 15 cm
DAG-503.12	Ø 2.5 mm X 15 cm
DAG-503.13	Ø 2.5 mm X 20 cm
DAG-503.14	Ø 3.0 mm X 15 cm
DAG-503.15	Ø 3.0 mm X 20 cm
DAG-503.16	Ø 3.0 mm X 25 cm

Luer lock	
DAG-504.11	Ø 2.0 mm X 15 cm
DAG-504.12	Ø 2.5 mm X 15 cm
DAG-504.13	Ø 2.5 mm X 20 cm
DAG-504.14	Ø 3.0 mm X 15 cm
DAG-504.15	Ø 3.0 mm X 20 cm
DAG-504.16	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)	
DAG-505.11	Ø 2.0 mm X 15 cm
DAG-505.12	Ø 2.5 mm X 15 cm
DAG-505.13	Ø 2.5 mm X 20 cm
DAG-505.14	Ø 3.0 mm X 15 cm
DAG-505.15	Ø 3.0 mm X 20 cm
DAG-505.16	Ø 3.0 mm X 25 cm

Power Cannula	
DAG-506.51	Ø 2 mm X 15 cm
DAG-506.52	Ø 2 mm X 20 cm
DAG-506.53	Ø 2 mm X 25 cm
DAG-506.54	Ø 2 mm X 30 cm
DAG-506.55	Ø 2 mm X 35 cm
DAG-506.56	Ø 3 mm X 15 cm
DAG-506.57	Ø 3 mm X 20 cm
DAG-506.58	Ø 3 mm X 25 cm
DAG-506.59	Ø 3 mm X 30 cm
DAG-506.60	Ø 3 mm X 35 cm
DAG-506.61	Ø 4 mm X 15 cm
DAG-506.62	Ø 4 mm X 20 cm
DAG-506.63	Ø 4 mm X 25 cm
DAG-506.64	Ø 4 mm X 30 cm
DAG-506.65	Ø 4 mm X 35 cm
DAG-506.66	Ø 5 mm X 15 cm
DAG-506.67	Ø 5 mm X 20 cm
DAG-506.68	Ø 5 mm X 25 cm
DAG-506.69	Ø 5 mm X 30 cm
DAG-506.70	Ø 5 mm X 35 cm
DAG-506.71	Ø 6 mm X 15 cm
DAG-506.72	Ø 6 mm X 20 cm
DAG-506.73	Ø 6 mm X 25 cm
DAG-506.74	Ø 6 mm X 30 cm
DAG-506.75	Ø 6 mm X 35 cm

Speciality Cannulas



Field Flap

With its elevated leading edge, this cannula is designed to open and expose soft tissue for more precise suctioning along the mandible and sub-mental areas.

60ccTommeyhub

DAG-501.21	Ø 2.5 mm X 15 cm
DAG-501.22	Ø 2.5 mm X 20 cm
DAG-501.23	Ø 3.0 mm X 15 cm
DAG-501.24	Ø 3.0 mm X 20 cm
DAG-501.25	Ø 3.0 mm X 25 cm
DAG-501.26	Ø 4.0 mm X 25 cm
DAG-501.27	Ø 4.0 mm X 30 cm
DAG-501.28	Ø 4.0 mm X 35 cm
DAG-501.29	Ø 4.5 mm X 30 cm
DAG-501.30	Ø 4.5 mm X 35 cm

Treaded Hub

DAG-502.21	Ø 3.0 mm X 15 cm
DAG-502.22	Ø 3.0 mm X 25 cm
DAG-502.23	Ø 4.0 mm X 15 cm
DAG-502.24	Ø 4.0 mm X 25 cm

Ultra Luer lock

DAG-503.21	Ø 2.0 mm X 15 cm
DAG-503.22	Ø 2.5 mm X 15 cm
DAG-503.23	Ø 2.5 mm X 20 cm
DAG-503.24	Ø 3.0 mm X 15 cm
DAG-503.25	Ø 3.0 mm X 20 cm
DAG-503.26	Ø 3.0 mm X 25 cm

Luer lock

DAG-504.21	Ø 2.0 mm X 15 cm
DAG-504.22	Ø 2.5 mm X 15 cm
DAG-504.23	Ø 2.5 mm X 20 cm
DAG-504.24	Ø 3.0 mm X 15 cm
DAG-504.25	Ø 3.0 mm X 20 cm
DAG-504.26	Ø 3.0 mm X 25 cm


Luer lock (Stainless Steel)

DAG-505.21	Ø 2.0 mm X 15 cm
DAG-505.22	Ø 2.5 mm X 15 cm
DAG-505.23	Ø 2.5 mm X 20 cm
DAG-505.24	Ø 3.0 mm X 15 cm
DAG-505.25	Ø 3.0 mm X 20 cm
DAG-505.26	Ø 3.0 mm X 25 cm

Power Cannula

DAG-506.101	Ø 2 mm X 15 cm
DAG-506.102	Ø 2 mm X 20 cm
DAG-506.103	Ø 2 mm X 25 cm
DAG-506.104	Ø 2 mm X 30 cm
DAG-506.105	Ø 2 mm X 35 cm
DAG-506.106	Ø 3 mm X 15 cm
DAG-506.107	Ø 3 mm X 20 cm
DAG-506.108	Ø 3 mm X 25 cm
DAG-506.109	Ø 3 mm X 30 cm
DAG-506.110	Ø 3 mm X 35 cm
DAG-506.111	Ø 4 mm X 15 cm
DAG-506.112	Ø 4 mm X 20 cm
DAG-506.113	Ø 4 mm X 25 cm
DAG-506.114	Ø 4 mm X 30 cm
DAG-506.115	Ø 4 mm X 35 cm
DAG-506.116	Ø 5 mm X 15 cm
DAG-506.117	Ø 5 mm X 20 cm
DAG-506.118	Ø 5 mm X 25 cm
DAG-506.119	Ø 5 mm X 30 cm
DAG-506.120	Ø 5 mm X 35 cm
DAG-506.121	Ø 6 mm X 15 cm
DAG-506.122	Ø 6 mm X 20 cm
DAG-506.123	Ø 6 mm X 25 cm
DAG-506.124	Ø 6 mm X 30 cm
DAG-506.125	Ø 6 mm X 35 cm

Speciality Cannulas



Ultra Luer Lock

The Super Luer Lock is on our Single-use instruments.

TumescentInfiltrator

This cannula is used specifically for the infusion of tumescent solution and for hydraulic dissection. the configuration of the ports allows even distribution of the fluid.

SL-601.01	Ø 2.0 mm X 12cm
SL-601.02	Ø 2.0 mm X 20cm

TonnardHarvester

A multiport harvester with 1mm elevated sharp microports. The 20 port configuration allows for harvesting both on the forward and backward strokes.

SL-602.01	Ø 3.0 mm X 20cm
SL-602.02	Ø 2.4 mm X 20cm

TonnardHarvester

Excellent for harvesting small aliquots of tissues.

SL-603.01	Ø 2.5 mm X 15cm
SL-603.02	Ø 2.0 mm X 15cm

CarrawayHarvester

Spiral three port design effective for multi-plane tissue harvesting.

SL-604.01	Ø 2.5 mm X 15cm
SL-604.02	Ø 2.0 mm X 10cm
SL-604.03	Ø 2.0 mm X 15cm
SL-604.04	Ø 1.5 mm X 8cm
SL-604.05	Ø 1.5 mm X 10cm

Injector

Round blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp injectors.

SL-605.01	Ø 1.0 mm X 5cm
SL-605.02	Ø 1.0 mm X 7cm
SL-605.03	Ø 1.0 mm X 5cm
SL-605.04	Ø 1.5 mm X 7cm
SL-605.05	Ø 1.0 mm X 4cm

SpoonTipInjector

Blunt injector designed with a more open port for rapid injection of larger volumes - a torso and buttock workhorse.

SL-606.01	Ø 3.0 mm X 20cm
SL-606.02	Ø 2.0 mm X 15cm

MillerSpeedHarvester

SL-607.01	Ø 2.0 mm X 15cm
-----------	-----------------

HOT
PRODUCT


Favourite

LIPOSUCTION ACCESSORIES

- AdaptorHandle
- LuerLockHandle
- NanoFilter
- AerobicTransfer
- AutologousNanoFatTransferLuertoLuerSetforLiposuction
- DecantStandAnaerobicTransfers
- SyringeHoldingStand
- JohnnieLok
- Crown
- PowerAdaptor
- PowerInjector
- FatSyringe
- Snaplok
- Two-WayValve
- 60ccHoldingSyringe
- IVSpike
- NeedleConnector


5253

Accessories *Liposuction*



AdaptorHandle

DAG-650.01



LuerLockHandle

DAG-650.02



ThreadedHandle

DAG-650.03



ThreadedHandle

DAG-650.04

Accessories *Liposuction*



NanoFilter

DAG-651.01 Introducing the patent-pending NanoTransfer. This proprietary closed system helps size adipose tissue so that it is injectable with 27g and 30g needles. NanoTransfer - taking fat transfer to the smallest level.



AnaerobicTransfers

Luer to Luer
DAG-652.12 1.2 mm
DAG-652.14 1.4 mm
DAG-652.24 2.4mm



AutologousNanoFatTransfer Luer-to-Luer Set-for-Liposuction

DAG-653.05 0.9 x 5
DAG-653.07 0.7 x 9
DAG-653.12 0.5 x 12

Accessories *Liposuction*



Decant Stand Anaerobic Transfers

Anaerobic transfers are used to prepare and transfer extracted tissue for re-injection. The Decant Stand design also allows for a natural gravity decant of extracted tissue.

DAG-654.01	60cc Toomey to 60cc Toomey
DAG-654.02	60cc Toomey to Luer Lock
DAG-654.03	60cc BD to Luer Lock
DAG-654.04	Luer Lock to Luer Lock



Syringe Holding Stands

DAG-655.01

Accessories *Liposuction*



SyringeStopper


DAG-656.01



JohnnieLok

The JohnnieLok holds the plunger under vacuum for easy syringe harvesting. The plunger is drawn back to the desired pressure level then twisted and locked securely.

DAG-657.01



Crown

A tiny device with a huge impact - Crowns change the centrifuge game. The Crown is designed to keep syringe contents in a virtually anaerobic atmosphere during centrifuging. Simply attach a Crown to the syringe plunger, harvest, then twist the plunger off. The Crown also has a bonus, self-sealing extractor port for cleanly removing contents from any layer of centrifuged tissue.

DAG-658.01	60cc Monoject
DAG-658.02	60cc BD
DAG-658.03	20cc BD
DAG-658.04	10cc BD

Accessories *Liposuction*

PowerAdaptor

The Tulip Power Adaptor modifies the syringe barrel for use with a liposuction aspirating machine.

DAG-659.01



PowerInjector

The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or 60cc Monoject syringe.

DAG-660.20 20cc
DAG-660.10 10cc
DAG-660.00 Universal
(for both 20 and 10cc
syringes)


Accessories *Liposuction*



FatSyringe

The fat syringe is an ideal tool to help surgeons for manual liposuction by creating negative pressure.

DAG-661.60 50/60 cc
DAG-661.20 20 cc
DAG-661.10 10 cc
DAG-661.05 5 cc
DAG-661.03 3 cc



SnapLoks

These syringe snap locks fit into the lunger flanges of the syringe. When the plunger is pulled out the lock snaps open on the barrel lip thus holding the syringe under vacuum.

DAG-662.60 60cc
DAG-662.50 50cc
DAG-662.35 35cc
DAG-662.20 20cc
DAG-662.10 10cc



Trocars

DAG-
669.01

Accessories *Liposuction*



Two-Way Valve.

Transfer the fat to the collection bag and take the fat out after treatment for injection.
DAG-663.01



3Way Stopcock

DAG-663.02



60cc Holding Stands

DAG-664.01



Cannula Adaptor

DAG-670.01

Accessories *Liposuction*



IV Spike

IV connector for use with reusable tubing.


DAG-665.01 stainless steel



Needle Connector

Needle connector for use with reusable tubing.

DAG-666.01 Length 2.3"



Skin Part

Skin Part for use in skin for pipe inserts.

DAG-667.01 6mm
DAG-667.02 5mm
DAG-667.03 4mm

Accessories *Liposuction*



Liposuctionjar

DAG-668.01



LIPOSUCTION *Workstation*